




Comparative Formulation

LB (Miller)

|  Powdered Media | LBP01 | LBP02 |
|-------------------------------------------------------------------------------------------------|------------|------------|
| | mg/L | mg/L |
| Agar, Bacteriological Type I | | 15000.0000 |
| Casein Digest Peptone | 10000.0000 | 10000.0000 |
| Sodium Chloride | 10000.0000 | 10000.0000 |
| Yeast Extract | 5000.0000 | 5000.0000 |



Liquid Media

LBL01

mg/L

Casein Digest Peptone

10000.0000

Sodium Chloride

10000.0000

Yeast Extract

5000.0000