



Liquid Media

TRL01

TRL02

TRL03

TRL04

mg/L

mg/L

mg/L

mg/L

D-Glucose, Anhydrous

1000.0000

1000.0000

EDTA, Disodium, Dihydrate

372.2000

372.2000

Phenol Red, Sodium Salt

17.0000

17.0000

Potassium Chloride

400.0000

400.0000

Potassium Phosphate, Monobasic, Anhydrous

60.0000

60.0000

Sodium Bicarbonate

350.0000

350.0000

Sodium Chloride

8000.0000

8500.0000

8500.0000

8000.0000

Sodium Phosphate Dibasic, Anhydrous

47.5000

47.5000

Tetrasodium-EDTA, Anhydrous

2205.8601

Trypsin Powder from Porcine Pancreas

2500.0000

5000.0000

25000.0000

500.0000