

Liquid Media	TRL01	TRL02	TRL03	TRL04	TRL05	TRL06
	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
D-Glucose, Anhydrous	1000.0000			1000.0000		
EDTA, Disodium, Dihydrate	372.2000			372.2000		
Phenol Red, Sodium Salt	17.0000			17.0000		
Potassium Chloride	400.0000			400.0000		
Potassium Phosphate, Monobasic, Anhydrous	60.0000			60.0000		
Sodium Bicarbonate	350.0000			350.0000		
Sodium Chloride	8000.0000	8500.0000	8500.0000	8000.0000		
Sodium Phosphate Dibasic, Anhydrous	47.5000			47.5000		
Tetrasodium-EDTA, Anhydrous		2205.8601				
Trypsin Powder from Porcine Pancreas	2500.0000	5000.0000	25000.0000	500.0000		