



# Formulations & Specifications

## PBS (10X), pH 7.2

**Product Number PBL07-6X1000ML**

Phosphate-buffered saline, 10X, pH 7.2. Store at room temperature.

Components	mg/L
Potassium Phosphate, Monobasic, Anhydrous (KH <sub>2</sub> PO <sub>4</sub> )	2100.0000
Sodium Chloride (NaCl)	90000.0000
Sodium Phosphate, Dibasic, Heptahydrate (Na <sub>2</sub> HPO <sub>4</sub> · 7H <sub>2</sub> O)	7260.0000