



Formulations & Specifications

PBS (10X), pH 7.2

Product Number PBL07-6X1000ML

Phosphate-buffered saline, 10X, pH 7.2. Store at room temperature.

| Components | mg/L |
|--|------------|
| Potassium Phosphate, Monobasic, Anhydrous (KH ₂ PO ₄) | 2100.0000 |
| Sodium Chloride (NaCl) | 90000.0000 |
| Sodium Phosphate, Dibasic, Heptahydrate (Na ₂ HPO ₄ · 7H ₂ O) | 7260.0000 |