

Gamborg's Vitamin Mixture (Powder)

Product Number GVP01-0.1LT

Contains the vitamins as described by Gamborg et al, 1968. Each liter equivalent of powder makes one liter of a 1000X vitamin stock solution. Store at 2° to 8°C.

| Components | mg/L |
|---|------------|
| Myo-Inositol (C ₆ H ₁₂ O ₆) | 100000.000 |
| | 0 |
| Nicotinic Acid (C ₆ H ₅ NO ₂) | 1000.0000 |
| Pyridoxine, Hydrochloride (C ₈ H ₁₁ NO ₃ . HCl) | 1000.0000 |
| Thiamine, Hydrochloride (C ₁₂ H ₁₇ CIN ₄ OS . HCl) | 10000.0000 |