

Lloyd & McCown's Woody Plant Basal Salts

Product Number WPP03-100LT

With macronutrients and micronutrients as described by
Lloyd and McCown, 1981. Store at 2° to 8°C.

Components	mg/L
Ammonium Nitrate (NH ₄ NO ₃)	400.0000
Boric Acid (H ₃ BO ₃)	6.2000
Calcium Chloride, Anhydrous (CaCl ₂)	72.4700
Calcium Nitrate, Tetrahydrate (Ca(NO ₃) ₂ · 4H ₂ O)	556.0000
Cupric Sulfate, Anhydrous (CuSO ₄)	0.1600
EDTA, Disodium, Dihydrate (C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ · 2H ₂ O)	37.3000
Ferrous Sulfate, Heptahydrate (FeSO ₄ ·7H ₂ O)	27.8500
Magnesium Sulfate, Anhydrous (MgSO ₄)	180.7000
Manganese Sulfate, Monohydrate (MnSO ₄ · H ₂ O)	22.3000
Molybdic Acid Sodium Salt, Dihydrate (Na ₂ MoO ₄ · 2H ₂ O)	0.2500
Potassium Phosphate, Monobasic, Anhydrous (KH ₂ PO ₄)	170.0000
Potassium Sulfate, Anhydrous (K ₂ SO ₄)	990.0000
Zinc Sulfate, Heptahydrate (ZnSO ₄ · 7H ₂ O)	8.6000