

## Lloyd & McCown's Woody Plant Basal Salts

**Product Number WPP03-100LT**

With macronutrients and micronutrients as described by  
Lloyd and McCown, 1981. Store at 2° to 8°C.

Components	mg/L
Ammonium Nitrate (NH <sub>4</sub> NO <sub>3</sub> )	400.0000
Boric Acid (H <sub>3</sub> BO <sub>3</sub> )	6.2000
Calcium Chloride, Anhydrous (CaCl <sub>2</sub> )	72.4700
Calcium Nitrate, Tetrahydrate (Ca(NO <sub>3</sub> ) <sub>2</sub> · 4H <sub>2</sub> O)	556.0000
Cupric Sulfate, Anhydrous (CuSO <sub>4</sub> )	0.1600
EDTA, Disodium, Dihydrate (C <sub>10</sub> H <sub>14</sub> N <sub>2</sub> Na <sub>2</sub> O <sub>8</sub> · 2H <sub>2</sub> O)	37.3000
Magnesium Sulfate, Anhydrous (MgSO <sub>4</sub> )	180.7000
Manganese Sulfate, Monohydrate (MnSO <sub>4</sub> · H <sub>2</sub> O)	22.3000
Molybdic Acid Sodium Salt, Dihydrate (Na <sub>2</sub> MoO <sub>4</sub> · 2H <sub>2</sub> O)	0.2500
Potassium Phosphate, Monobasic, Anhydrous (KH <sub>2</sub> PO <sub>4</sub> )	170.0000
Potassium Sulfate, Anhydrous (K <sub>2</sub> SO <sub>4</sub> )	990.0000
Zinc Sulfate, Heptahydrate (ZnSO <sub>4</sub> · 7H <sub>2</sub> O)	8.6000