

Lloyd & McCown's Woody Plant Basal Salts

Product Number WPP03-100LT

With macronutrients and micronutrients as described by
Lloyd and McCown, 1981. Store at 2° to 8°C.

| Components | mg/L |
|---|----------|
| Ammonium Nitrate (NH ₄ NO ₃) | 400.0000 |
| Boric Acid (H ₃ BO ₃) | 6.2000 |
| Calcium Chloride, Anhydrous (CaCl ₂) | 72.4700 |
| Calcium Nitrate, Tetrahydrate (Ca(NO ₃) ₂ · 4H ₂ O) | 556.0000 |
| Cupric Sulfate, Anhydrous (CuSO ₄) | 0.1600 |
| EDTA, Disodium, Dihydrate (C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ · 2H ₂ O) | 37.3000 |
| Magnesium Sulfate, Anhydrous (MgSO ₄) | 180.7000 |
| Manganese Sulfate, Monohydrate (MnSO ₄ · H ₂ O) | 22.3000 |
| Molybdic Acid Sodium Salt, Dihydrate (Na ₂ MoO ₄ · 2H ₂ O) | 0.2500 |
| Potassium Phosphate, Monobasic, Anhydrous (KH ₂ PO ₄) | 170.0000 |
| Potassium Sulfate, Anhydrous (K ₂ SO ₄) | 990.0000 |
| Zinc Sulfate, Heptahydrate (ZnSO ₄ · 7H ₂ O) | 8.6000 |