



Product Catalog **2009**

Plant Tissue Culture
Antibiotics/Antimycotics

Product Catalog

Plant Tissue Culture Antibiotics / Antimycotics (Plant)

Amphotericin B

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Amphotericin B Solution/ Contains 250µg/ml amphotericin B and 250µg/ml sodium deoxycholate. Effective against fungi and yeasts. Sterile solution.	Liquid	100	ML	ABL01-100ML	Inquire	Inquire

Antibiotic-Antimycotic

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Antibiotic-Antimycotic Solution (100X)/ Contains 10,000 units of penicillin, 10,000µg streptomycin, and 25µg amphotericin B in a 0.85% saline solution. Effective against bacteria, fungi, and yeasts. Sterile solution.	Liquid	100	ML	ABL02-100ML	Inquire	Inquire

Carbenicillin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Carbenicillin, Disodium Salt/ >=90%	Powder	10	G	C033-10GM	Inquire	Inquire
Carbenicillin, Disodium Salt/ >=90%	Powder	25	G	C033-25GM	Inquire	Inquire

Cefotaxime

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Cefotaxime, Sodium Salt/ >=95%	Powder	10	G	C032-10GM	Inquire	Inquire
Cefotaxime, Sodium Salt/ >=95%	Powder	25	G	C032-25GM	Inquire	Inquire

Product Catalog

Plant Tissue Culture Antibiotics / Antimycotics (Plant)

Gentamicin						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Gentamicin Solution (50mg/ml)/ Contains 50mg/ml gentamicin. Sterile solution.	Liquid	10	ML	ABL03-10ML	Inquire	Inquire
Gentamicin Sulfate/ >=590µg/mg	Powder	5	G	G020-5GM	Inquire	Inquire
Gentamicin Sulfate/ >=590µg/mg	Powder	10	G	G020-10GM	Inquire	Inquire
Gentamicin Sulfate/ >=590µg/mg	Powder	25	G	G020-25GM	Inquire	Inquire

Hygromycin B						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Hygromycin B/ >=60% (HPAE)	Powder	1	G	H010-1GM	Inquire	Inquire
Hygromycin B/ >=60% (HPAE)	Powder	5	G	H010-5GM	Inquire	Inquire

Kanamycin						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Kanamycin Sulfate/ >=670 units/mg	Powder	10	G	K003-10GM	Inquire	Inquire
Kanamycin Sulfate/ >=670 units/mg	Powder	25	G	K003-25GM	Inquire	Inquire

Penicillin G						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Penicillin G, Sodium Salt/ >=96%; >=590µg/mg	Powder	25	G	P056-25GM	Inquire	Inquire
Penicillin G, Sodium Salt/ >=96%; >=590µg/mg	Powder	100	G	P056-100GM	Inquire	Inquire
Penicillin G, Sodium Salt/ >=96%; >=590µg/mg	Powder	500	G	P056-500GM	Inquire	Inquire

Product Catalog

Plant Tissue Culture Antibiotics / Antimycotics (Plant)

Penicillin-Streptomycin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Penicillin-Streptomycin Solution (100X)/ Contains 10,000 units of penicillin and 10,000µg streptomycin in a 0.85% saline solution. Sterile solution.	Liquid	100	ML	PSL01-100ML	Inquire	Inquire

Plant Preservative Mixture (PPM)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Plant Preservative Mixture (PPM)/ Broad-spectrum biocide/fungicide for plant tissue culture.	Liquid	10	ML	PPL33-10ML	Inquire	Inquire
Plant Preservative Mixture (PPM)/ Broad-spectrum biocide/fungicide for plant tissue culture.	Liquid	30	ML	PPL33-30ML	Inquire	Inquire
Plant Preservative Mixture (PPM)/ Broad-spectrum biocide/fungicide for plant tissue culture.	Liquid	100	ML	PPL33-100ML	Inquire	Inquire
Plant Preservative Mixture (PPM)/ Broad-spectrum biocide/fungicide for plant tissue culture.	Liquid	500	ML	PPL33-500ML	Inquire	Inquire

Streptomycin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Streptomycin Sulfate/ >=650 units/mg	Powder	25	G	S041-25GM	Inquire	Inquire
Streptomycin Sulfate/ >=650 units/mg	Powder	100	G	S041-100GM	Inquire	Inquire
Streptomycin Sulfate/ >=650 units/mg	Powder	500	G	S041-500GM	Inquire	Inquire

Product Catalog
Plant Tissue Culture Antibiotics / Antimycotics (Plant)

Timentin						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Timentin/ 15:1 Ticarcillin Disodium Salt / Potassium Clavulanate Mixture	Powder	1	G	T034-1GM	Inquire	Inquire
Timentin/ 15:1 Ticarcillin Disodium Salt / Potassium Clavulanate Mixture	Powder	5	G	T034-5GM	Inquire	Inquire
Timentin/ 15:1 Ticarcillin Disodium Salt / Potassium Clavulanate Mixture	Powder	25	G	T034-25GM	Inquire	Inquire

Plant Tissue Culture
Biochemicals and Other Reagents

Product Catalog

Plant Tissue Culture Biochemicals and Reagents

Ammonium Nitrate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Ammonium Nitrate/ CAS# 6484-52-2; >=98%; HAZARDOUS MATERIAL	Powder	1	KG	A015-1KG	Inquire	Inquire
Ammonium Nitrate/ CAS# 6484-52-2; >=98%; HAZARDOUS MATERIAL	Powder	5	KG	A015-5KG	Inquire	Inquire
Ammonium Nitrate/ CAS# 6484-52-2; >=98%; HAZARDOUS MATERIAL	Powder	10	KG	A015-10KG	Inquire	Inquire

Boric Acid

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Boric Acid/ CAS# 10043-35-3; >=99.5%;	Powder	500	G	B003-500GM	Inquire	Inquire
Boric Acid/ CAS# 10043-35-3; >=99.5%;	Powder	2.5	KG	B003-2.5KG	Inquire	Inquire

Calcium Chloride

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Calcium Chloride, Anhydrous/ CAS# 10043-52-4; >=96%	Powder	50	G	C021-500GM	Inquire	Inquire
Calcium Chloride, Anhydrous/ CAS# 10043-52-4; >=96%	Powder	2.5	KG	C021-2.5KG	Inquire	Inquire

Calcium Nitrate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Calcium Nitrate, 4H2O/ CAS# 13477-34-4; >=99%; HAZARDOUS MATERIAL	Powder	1	KG	C017-1KG	Inquire	Inquire
Calcium Nitrate, 4H2O/ CAS# 13477-34-4; >=99%; HAZARDOUS MATERIAL	Powder	5	KG	C017-5KG	Inquire	Inquire

Product Catalog

Plant Tissue Culture Biochemicals and Reagents

Charcoal

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Charcoal, Activated & Acid Washed/ CAS# 7440-44-0	Powder	1	KG	C006-1KG	Inquire	Inquire

Choline Chloride

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Choline Chloride/ CAS# 67-48-1; >=98%	Powder	100	G	C013-100GM	Inquire	Inquire
Choline Chloride/ CAS# 67-48-1; >=98%	Powder	500	G	C013-500GM	Inquire	Inquire

Cobalt Chloride

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Cobalt Chloride, 6H2O/ CAS# 7791-13-1; >=98%	Powder	5	G	C007-5GM	Inquire	Inquire
Cobalt Chloride, 6H2O/ CAS# 7791-13-1; >=98%	Powder	25	G	C007-25GM	Inquire	Inquire

Cupric Sulfate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Cupric Sulfate, 5H2O/ CAS# 7758-99-8; >=98%	Powder	25	G	C009-25GM	Inquire	Inquire
Cupric Sulfate, 5H2O/ CAS# 7758-99-8; >=98%	Powder	500	G	C009-500GM	Inquire	Inquire

D-Glucose

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
D-Glucose, Anhydrous/ CAS# 50-99-7	Powder	500	G	G013-500GM	Inquire	Inquire
D-Glucose, Anhydrous/ CAS# 50-99-7	Powder	2.5	KG	G013-2.5KG	Inquire	Inquire
D-Glucose, Anhydrous/ CAS# 50-99-7	Powder	10	KG	G013-10KG	Inquire	Inquire

Dichloroisocyanuric Acid

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Dichloroisocyanuric Acid, Sodium Salt (NaDCC)/ CAS: 2893-78-9 EINECCS: 220-767-7 >=96% HAZARDOUS MATERIAL	Powder	1	KG	D011-1KG	Inquire	Inquire
Dichloroisocyanuric Acid, Sodium Salt/ CAS: 2893-78-9 EINECCS: 220-767-7 >=96% HAZARDOUS MATERIAL	Powder	500	GM	D011-500GM	Inquire	Inquire

EDTA

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Na₂ EDTA, 2H₂O/ CAS# 6381-92-6; >=99%	Powder	100	G	E004-100GM	Inquire	Inquire
Na₂ EDTA, 2H₂O/ CAS# 6381-92-6; >=99%	Powder	500	G	E004-500GM	Inquire	Inquire
Na₂ EDTA, 2H₂O/ CAS# 6381-92-6; >=99%	Powder	2.5	KG	E004-2.5KG	Inquire	Inquire

Ferrous Sulfate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Ferrous Sulfate, 7H₂O/ CAS# 7782-63-0; >=99.5%	Powder	100	G	F006-100GM	Inquire	Inquire
Ferrous Sulfate, 7H₂O/ CAS# 7782-63-0; >=99.5%	Powder	500	G	F006-500GM	Inquire	Inquire
Ferrous Sulfate, 7H₂O/ CAS# 7782-63-0; >=99.5%	Powder	2.5	KG	F006-2.5KG	Inquire	Inquire

Product Catalog

Plant Tissue Culture Biochemicals and Reagents

Glycine

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Glycine/ CAS# 56-40-6; >=98.5%	Powder	100	G	G007-100GM	Inquire	Inquire
Glycine/ CAS# 56-40-6; >=98.5%	Powder	1	KG	G007-1KG	Inquire	Inquire

Magnesium Sulfate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Magnesium Sulfate, 7H2O/ CAS# 10034-99-8; >=98%	Powder	500	G	M020-500GM	Inquire	Inquire
Magnesium Sulfate, 7H2O/ CAS# 10034-99-8; >=98%	Powder	1	KG	M020-1KG	Inquire	Inquire
Magnesium Sulfate, 7H2O/ CAS# 10034-99-8; >=98%	Powder	2.5	KG	M020-2.5KG	Inquire	Inquire
Magnesium Sulfate, 7H2O/ CAS# 10034-99-8; >=98%	Powder	10	KG	M020-10KG	Inquire	Inquire
Magnesium Sulfate, Anhydrous/ CAS# 7487-88-9; >=98%	Powder	500	G	M022-500GM	Inquire	Inquire
Magnesium Sulfate, Anhydrous/ CAS# 7487-88-9; >=98%	Powder	2.5	KG	M022-2.5KG	Inquire	Inquire
Magnesium Sulfate, Anhydrous/ CAS# 7487-88-9; >=98%	Powder	10	KG	M022-10KG	Inquire	Inquire

Manganese Sulfate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Manganese Sulfate, H2O/ CAS# 10034-96-5; >=98%	Powder	500	G	M006-500GM	Inquire	Inquire
Manganese Sulfate, H2O/ CAS# 10034-96-5; >=98%	Powder	2.5	KG	M006-2.5KG	Inquire	Inquire

Product Catalog

Plant Tissue Culture Biochemicals and Reagents

MES Buffer

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
MES Buffer/ CAS# 4432-31-9; >=99%	Powder	100	G	M009-100GM	Inquire	Inquire
MES Buffer/ CAS# 4432-31-9; >=99%	Powder	500	G	M009-500GM	Inquire	Inquire

Myo-Inositol

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
myo-Inositol/ CAS# 87-89-8; >=99%	Powder	500	G	M014-500GM	Inquire	Inquire
myo-Inositol/ CAS# 87-89-8; >=99%	Powder	2.5	KG	M014-2.5KG	Inquire	Inquire
myo-Inositol/ CAS# 87-89-8; >=99%	Powder	10	KG	M014-10KG	Inquire	Inquire

Nicotinic Acid

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Nicotinic Acid/ CAS# 59-67-6; >=98%	Powder	100	G	N004-100GM	Inquire	Inquire
Nicotinic Acid/ CAS# 59-67-6; >=98%	Powder	500	G	N004-500GM	Inquire	Inquire

Potassium Iodide

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Potassium Iodide/ CAS# 7681-11-0; >=99%	Powder	100	G	P011-100GM	Inquire	Inquire
Potassium Iodide/ CAS# 7681-11-0; >=99%	Powder	500	G	P011-500GM	Inquire	Inquire

Product Catalog

Plant Tissue Culture Biochemicals and Reagents

Potassium Nitrate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Potassium Nitrate/ CAS# 7757-79-1; >=99%; HAZARDOUS MATERIAL	Powder	1	KG	P012-1KG	Inquire	Inquire
Potassium Nitrate/ CAS# 7757-79-1; >=99%; HAZARDOUS MATERIAL	Powder	5	KG	P012-2.5KG	Inquire	Inquire
Potassium Nitrate/ CAS# 7757-79-1; >=99%; HAZARDOUS MATERIAL	Powder	10	KG	P012-10KG	Inquire	Inquire

Potassium Phosphate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Potassium Phosphate Dibasic, Anhydrous/ CAS# 7758-11-4; >=98%	Powder	500	G	P021-500GM	Inquire	Inquire
Potassium Phosphate Dibasic, Anhydrous/ CAS# 7758-11-4; >=98%	Powder	2.5	KG	P021-2.5KG	Inquire	Inquire
Potassium Phosphate Monobasic, Anhydrous/ CAS# 7778-77-0; >=99%	Powder	500	G	P031-500GM	Inquire	Inquire
Potassium Phosphate Monobasic, Anhydrous/ CAS# 7778-77-0; >=99%	Powder	2.5	KG	P031-2.5KG	Inquire	Inquire

Pyridoxine

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Pyridoxine HCl/ CAS# 58-56-0; >=98%	Powder	25	G	P018-25GM	Inquire	Inquire
Pyridoxine HCl/ CAS# 58-56-0; >=98%	Powder	25	G	P018-25GM	Inquire	Inquire

Product Catalog

Plant Tissue Culture Biochemicals and Reagents

Sodium Phosphate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Sodium Phosphate Dibasic/ CN #7558-79-4	Powder	500		S018-500GM	Inquire	Inquire

Sucrose

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Sucrose/ CAS# 57-50-1	Powder	500	G	S011-500GM	Inquire	Inquire
Sucrose/ CAS# 57-50-1	Powder	2.5	KG	S011-2.5KG	Inquire	Inquire
Sucrose/ CAS# 57-50-1	Powder	10	KG	S011-10KG	Inquire	Inquire

Thiamine

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Thiamine HCl/ CAS# 67-03-8; >=99%	Powder	25	G	T009-25GM	Inquire	Inquire
Thiamine HCl/ CAS# 67-03-8; >=99%	Powder	25	G	T009-25GM	Inquire	Inquire

Water, Cell Culture Grade

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Purified Water/ 18 mega-ohm, purified water. Cell culture tested. Sterile Filtered	Liquid	1000	ML	WAL01-1000ML	Inquire	Inquire

Zinc Sulfate

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Zinc Sulfate, 7H2O/ CAS# 7446-20-0; >=99%	Powder	100	G	Z006-100GM	Inquire	Inquire
Zinc Sulfate, 7H2O/ CAS# 7446-20-0; >=99%	Powder	500	G	Z006-500GM	Inquire	Inquire

Plant Tissue Culture
Bioreactors

Liquid Lab						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
<p>Liquid Lab Vent Patch/ Breathable patch with a 0.2 micron membrane that enables gas exchange in virtually any type of tissue culture vessel. The patch is bonded to a non-woven Polyester® for strength and is durable enough to withstand multiple passes through an autoclave. The Liquid Lab Patch® is made to withstand the rigors of a tissue culture lab with a Gurley airflow of 17 seconds/100cc/SI, water intrusion pressure of +50psi and a burst strength of 21psi. Strip of 10 patches.</p>		10	PC	SSVENT-10PC	Inquire	Inquire
<p>Liquid Lab Vessel/ 10.5"x 4"x 4"; Polycarbonate; Comes with microporous patch and a septaport</p>		1	PCS	SSV01-1PC	Inquire	Inquire
<p>M200/ Unit is designed to be used in large scale production labs. Standing eight feet tall, the M200 can hold up to two hundred vessels, on ten shelves. Unit comes standard with ten shelves, ten high output lights, control panel, and two hundred Southern Sun Liquid Lab™ vessels with either vent patch and septa-port or two vent-patches. Dimensions are 8' high x 5' long x 2.5' wide</p>		1	PCS	M200-1PC	Inquire	Inquire

Product Catalog
Plant Tissue Culture Bioreactors

<p>M34/ The M34 unit is designed for small scale production and research. The size and weight of the unit permits it to be placed on shelf or tabletop. This unit comes standard with two rocking shelves; two high output lights; control panel and thirty two Southern Sun Liquid Lab™ vessels either vent patch and septa-port or two vent-patches. Dimensions are 3.5' high x 4' long x 2.5' wide</p>		1	PCS	M34-1PC	Inquire	Inquire
---	--	---	-----	---------	---------	---------

Plant Tissue Culture
Gelling Agents

Agars						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Micropropagation Agar-Type I/ Suitable for most tissue culture applications; Use at 5 - 9g/L	Powder	500	G	A038-500GM	Inquire	Inquire
Micropropagation Agar-Type I/ Suitable for most tissue culture applications; Use at 5 - 9g/L	Powder	1	KG	A038-1KG	Inquire	Inquire
Micropropagation Agar-Type I/ Suitable for most tissue culture applications; Use at 5 - 9g/L	Powder	2.5	KG	A038-2.5KG	Inquire	Inquire
Micropropagation Agar-Type I/ Suitable for most tissue culture applications; Use at 5 - 9g/L	Powder	10	KG	A038-10KG	Inquire	Inquire
Micropropagation Agar-Type II/ Suitable for more sensitive applications; Use at 5 - 9g/L	Powder	500	G	A037-500GM	Inquire	Inquire
Micropropagation Agar-Type II/ Suitable for more sensitive applications; Use at 5 - 9g/L	Powder	1	KG	A037-1KG	Inquire	Inquire
Micropropagation Agar-Type II/ Suitable for more sensitive applications; Use at 5 - 9g/L	Powder	2.5	KG	A037-2.5KG	Inquire	Inquire
Micropropagation Agar-Type II/ Suitable for more sensitive applications; Use at 5 - 9g/L	Powder	10	KG	A037-10KG	Inquire	Inquire
Phytoblend/ Proprietary blend of agars (compare to Phytagar®) Use at 7 - 9g/L	Powder	500	G	PTP01-500GM	Inquire	Inquire
Phytoblend/ Proprietary blend of agars (compare to Phytagar®) Use at 7 - 9g/L	Powder	1	KG	PTP01-1KG	Inquire	Inquire

Product Catalog

Plant Tissue Culture Gelling Agents

Phytoblend/ Proprietary blend of agars (compare to Phytagar®) Use at 7 - 9g/L	Powder	2.5	KG	PTP01-2.5KG	Inquire	Inquire
Phytoblend/ Proprietary blend of agars (compare to Phytagar®) Use at 7 - 9g/L	Powder	10	KG	PTP01-10KG	Inquire	Inquire
UltaPure Agar/ Purified agar suitable for very sensitive applications; Use at 5.5 - 8.5g/L	Powder	100	G	A046-100GM	Inquire	Inquire
UltaPure Agar/ Purified agar suitable for very sensitive applications; Use at 5.5 - 8.5g/L	Powder	500	G	A046-500GM	Inquire	Inquire
UltaPure Agar/ Purified agar suitable for very sensitive applications; Use at 5.5 - 8.5g/L	Powder	1	KG	A046-1KG	Inquire	Inquire
UltaPure Agar/ Purified agar suitable for very sensitive applications; Use at 5.5 - 8.5g/L	Powder	2.5	KG	A046-2.5KG	Inquire	Inquire

Product Catalog

Plant Tissue Culture Gelling Agents

Other Gelling Agents						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Gellan Gum/ Produces a transparent gel; Use at 2 - 4g/L	Powder	500	G	G017-500GM	Inquire	Inquire
Gellan Gum/ Produces a transparent gel; Use at 2 - 4g/L	Powder	1	KG	G017-1KG	Inquire	Inquire
Gellan Gum/ Produces a transparent gel; Use at 2 - 4g/L	Powder	2.5	KG	G017-2.5KG	Inquire	Inquire
Gellan Gum/ Produces a transparent gel; Use at 2 - 4g/L	Powder	10	KG	G017-10KG	Inquire	Inquire
Gelzan/ Highest grade gellan gum available (Formerly named Gelrite®); Produces a very transparent gel; Use at 2 - 4g/L	Powder	500	G	G024-500GM	Inquire	Inquire
Gelzan/ Highest grade gellan gum available (Formerly named Gelrite®); Produces a very transparent gel; Use at 2 - 4g/L	Powder	1	KG	G024-1KG	Inquire	Inquire
Gelzan/ Highest grade gellan gum available (Formerly named Gelrite®); Produces a very transparent gel; Use at 2 - 4g/L	Powder	2.5	KG	G024-2.5KG	Inquire	Inquire
Gelzan/ Highest grade gellan gum available (Formerly named Gelrite®); Produces a very transparent gel; Use at 2 - 4g/L	Powder	10	KG	G024-10KG	Inquire	Inquire
TC Gel/ Proprietary blend of agar and Gelzan; Use at 4 - 6g/L	Powder	500	G	PTP02-500GM	Inquire	Inquire
TC Gel/ Proprietary blend of agar and Gelzan; Use at 4 - 6g/L	Powder	1	KG	PTP02-1KG	Inquire	Inquire
TC Gel/ Proprietary blend of agar and Gelzan; Use at 4 - 6g/L	Powder	2.5	KG	PTP02-2.5KG	Inquire	Inquire
TC Gel/ Proprietary blend of agar and Gelzan; Use at 4 - 6g/L	Powder	10	KG	PTP02-10KG	Inquire	Inquire

Plant Tissue Culture
Plant Growth Regulators

Product Catalog

Plant Tissue Culture Plant Growth Regulators

(+)-cis,trans-Abscisic Acid (ABA)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
(+)-cis,trans-Abscisic Acid (ABA)/ >98%	Powder	100	MG	A036-100MG	Inquire	Inquire
(+)-cis,trans-Abscisic Acid (ABA)/ >98%	Powder	250	MG	A036-250MG	Inquire	Inquire
(+)-cis,trans-Abscisic Acid (ABA)/ >98%	Powder	500	MG	A036-500MG	Inquire	Inquire

2-carboxyphenyl 3-phenylpropane 1,3-dione (CPD)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
2-carboxyphenyl 3-phenylpropane 1, 3-dione (CPD)/ Auxin inhibitor	Powder	100	MG	C043-100MG	Inquire	Inquire
2-carboxyphenyl 3-phenylpropane 1, 3-dione (CPD)/ Auxin inhibitor	Powder	500	MG	C043-500MG	Inquire	Inquire
2-carboxyphenyl 3-phenylpropane 1, 3-dione (CPD)/ Auxin inhibitor	Powder	1	G	C043-1GM	Inquire	Inquire

6-(g,g-dimethylallylamino purine) (2iP)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
6-(g,g-dimethylallylamino purine) (2iP)/ >99%	Powder	1	G	D006-1GM	Inquire	Inquire
6-(g,g-dimethylallylamino purine) (2iP)/ >99%	Powder	10	G	D006-10GM	Inquire	Inquire
6-(g,g-dimethylallylamino purine) (2iP)/ >99%	Powder	25	G	D006-25GM	Inquire	Inquire

Product Catalog

Plant Tissue Culture Plant Growth Regulators

6-Benzylaminopurine (6-BA, BAP)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
6-Benzylaminopurine (6-BA, BA, BAP)/ >99%	Powder	5	G	B001-5GM	Inquire	Inquire
6-Benzylaminopurine (6-BA, BA, BAP)/ >99%	Powder	25	G	B001-25GM	Inquire	Inquire
6-Benzylaminopurine (6-BA, BA, BAP)/ >99%	Powder	100	G	B001-100GM	Inquire	Inquire
6-Benzylaminopurine (6-BA, BA, BAP)/ >99%	Powder	1	KG	B001-1KG	Inquire	Inquire

6-Furfurylaminopurine (Kinetin)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
6-Furfurylaminopurine (Kinetin)/ >99%	Powder	1	G	K001-1GM	Inquire	Inquire
6-Furfurylaminopurine (Kinetin)/ >99%	Powder	5	G	K001-5GM	Inquire	Inquire
6-Furfurylaminopurine (Kinetin)/ >99%	Powder	25	G	K001-25GM	Inquire	Inquire
6-Furfurylaminopurine (Kinetin)/ >99%	Powder	100	G	K001-100GM	Inquire	Inquire

GA4+7

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
GA4+7/ >90%	Powder	1	G	G026-1GM	Inquire	Inquire
GA4+7/ >90%	Powder	5	G	G026-5GM	Inquire	Inquire

Product Catalog

Plant Tissue Culture Plant Growth Regulators

Gibberellic Acid (GA3)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Gibberellic Acid (GA3)/ >90%	Powder	1	G	G001-1GM	Inquire	Inquire
Gibberellic Acid (GA3)/ >90%	Powder	5	G	G001-5GM	Inquire	Inquire
Gibberellic Acid (GA3)/ >90%	Powder	25	G	G001-25GM	Inquire	Inquire
Gibberellic Acid (GA3)/ >90%	Powder	100	G	G001-100GM	Inquire	Inquire

Indole-3-acetic Acid (IAA)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Indole-3-acetic Acid (IAA)/ >99%	Powder	5	G	I001-5GM	Inquire	Inquire
Indole-3-acetic Acid (IAA)/ >99%	Powder	25	G	I001-25GM	Inquire	Inquire
Indole-3-acetic Acid (IAA)/ >99%	Powder	100	G	I001-100GM	Inquire	Inquire

Indole-3-butyric Acid (IBA)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Indole-3-butyric Acid (IBA)/ >99%	Powder	5	G	I002-5GM	Inquire	Inquire
Indole-3-butyric Acid (IBA)/ >99%	Powder	25	G	I002-25GM	Inquire	Inquire
Indole-3-butyric Acid (IBA)/ >99%	Powder	100	G	I002-100GM	Inquire	Inquire
Indole-3-butyric Acid (IBA)/ >99%	Powder	500	G	I002-500GM	Inquire	Inquire
Indole-3-butyric Acid (IBA)/ >99%	Powder	1	KG	I002-1KG	Inquire	Inquire

Product Catalog

Plant Tissue Culture Plant Growth Regulators

Indole-3-propionic Acid (IPA)						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Indole-3-propionic Acid (IPA)/ >99%	Powder	5	G	I009-5GM	Inquire	Inquire
Indole-3-propionic Acid (IPA)/ >99%	Powder	25	G	I009-25GM	Inquire	Inquire
Indole-3-propionic Acid (IPA)/ >99%	Powder	100	G	I009-100GM	Inquire	Inquire
Indole-3-propionic Acid (IPA)/ >99%	Powder	500	G	I009-500GM	Inquire	Inquire

meta-Topolin						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
meta-Topolin/ CAS# 75737-38-1; C12H11N5O; MW 241.25 Plant tissue culture tested	Powder	100	MG	T039-100MG	Inquire	Inquire
meta-Topolin/ CAS# 75737-38-1; C12H11N5O; MW 241.25 Plant tissue culture tested	Powder	500	MG	T039-500MG	Inquire	Inquire
meta-Topolin/ CAS# 75737-38-1; C12H11N5O; MW 241.25 Plant tissue culture tested	Powder	1	G	T039-1GM	Inquire	Inquire

Naphthaleneacetic Acid (NAA)						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Naphthaleneacetic Acid (NAA)/ >98%	Powder	25	G	N001-25GM	Inquire	Inquire
Naphthaleneacetic Acid (NAA)/ >98%	Powder	100	G	N001-100GM	Inquire	Inquire
Naphthaleneacetic Acid (NAA)/ >98%	Powder	1	KG	N001-1KG	Inquire	Inquire
Naphthaleneacetic Acid (NAA) Solution/ 1 mg/ml in H2O. Sterile Filtered	Liquid	100	ML	NSL01-100ML	Inquire	Inquire

Product Catalog

Plant Tissue Culture Plant Growth Regulators

Paclobutrazol

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Paclobutrazol/ >=95%	Powder	25	G	P054-25GM	Inquire	Inquire
Paclobutrazol/ >=95%	Powder	100	G	P054-100GM	Inquire	Inquire

Thidiazuron (TDZ)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Thidiazuron (TDZ)/ >98%	Powder	100	MG	T017-100MG	Inquire	Inquire
Thidiazuron (TDZ)/ >98%	Powder	100	MG	T017-100MG	Inquire	Inquire
Thidiazuron (TDZ)/ >98%	Powder	500	MG	T017-500MG	Inquire	Inquire
Thidiazuron (TDZ)/ >98%	Powder	1	G	T017-1G	Inquire	Inquire

Zeatin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
trans-Zeatin/ >99%	Powder	100	MG	Z007-100MG	Inquire	Inquire
trans-Zeatin/ >99%	Powder	500	MG	Z007-500MG	Inquire	Inquire
trans-Zeatin/ >99%	Powder	1	G	Z007-1G	Inquire	Inquire
trans-Zeatin Solution/ Contains 1 mg/ml in H ₂ O. Sterile filtered.	Liquid	10	ML	ZSL01-10ML	Inquire	Inquire

Product Catalog
Plant Tissue Culture Plant Growth Regulators

Zeatin Riboside						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
trans-Zeatin Riboside/ >99%	Powder	50	MG	Z008-50MG	Inquire	Inquire
trans-Zeatin Riboside/ >99%	Powder	100	MG	Z008-100MG	Inquire	Inquire
trans-Zeatin Riboside/ >99%	Powder	250	MG	Z008-250MG	Inquire	Inquire
trans-Zeatin Riboside/ >99%	Powder	500	MG	Z008-500MG	Inquire	Inquire

Plant Tissue Culture Media

Anderson's Rhododendron						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Anderson's Rhododendron/ With macro- and micronutrients.	Powder	1	L	ARP01-1LT	Inquire	Inquire
Anderson's Rhododendron/ With macro- and micronutrients.	Powder	10	L	ARP01-10LT	Inquire	Inquire
Anderson's Rhododendron/ With macro- and micronutrients.	Powder	50	L	ARP01-50LT	Inquire	Inquire
Anderson's Rhododendron with Vitamins/ With macro- and micronutrients and LS vitamins.	Powder	1	L	ARP05-1LT	Inquire	Inquire
Anderson's Rhododendron with Vitamins/ With macro- and micronutrients and LS vitamins.	Powder	10	L	ARP05-10LT	Inquire	Inquire
Anderson's Rhododendron with Vitamins/ With macro- and micronutrients and LS vitamins.	Powder	50	L	ARP05-50LT	Inquire	Inquire

Cereal Crops						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Cereal Crops Medium/ With macro- and micronutrients, amino acids, and carbohydrates.	Powder	1	L	CCP01-1LT	Inquire	Inquire
Cereal Crops Medium/ With macro- and micronutrients, amino acids, and carbohydrates.	Powder	5	L	CCP01-5LT	Inquire	Inquire
Cereal Crops Medium/ With macro- and micronutrients, amino acids, and carbohydrates.	Powder	10	L	CCP01-10LT	Inquire	Inquire

Chee & Pool						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Chee & Pool (C2D) Vitis Medium/ With macro- and micronutrients	Powder	1	L	CEP01-1LT	Inquire	Inquire
Chee & Pool (C2D) Vitis Medium/ With macro- and micronutrients	Powder	10	L	CEP01-10LT	Inquire	Inquire
Chee & Pool (C2D) Vitis Medium/ With macro- and micronutrients	Powder	50	L	CEP01-50LT	Inquire	Inquire
Chee & Pool (C2D) Vitis Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	1	L	CEP02-1LT	Inquire	Inquire
Chee & Pool (C2D) Vitis Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	50	L	CEP02-50LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Chu's						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Chu's N-6 Medium/ With macro- and micronutrients.	Powder	1	L	CHP01-1LT	Inquire	Inquire
Chu's N-6 Medium/ With macro- and micronutrients.	Powder	10	L	CHP01-10LT	Inquire	Inquire
Chu's N-6 Medium/ With macro- and micronutrients.	Powder	50	L	CHP01-50LT	Inquire	Inquire
Chu's N-6 Medium with Vitamins/ With macro- and micronutrients and vitamins.	Powder	1	L	CHP03-1LT	Inquire	Inquire
Chu's N-6 Medium with Vitamins/ With macro- and micronutrients and vitamins.	Powder	10	L	CHP03-10LT	Inquire	Inquire
Chu's N-6 Medium with Vitamins/ With macro- and micronutrients and vitamins.	Powder	50	L	CHP03-50LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

CLC/Ipomoea						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Callus Proliferation Medium/ With macro- and micronutrients and vitamins.	Powder	0.5	L	CPP01-0.5LT	Inquire	Inquire
Callus Proliferation Medium/ With macro- and micronutrients and vitamins.	Powder	1	L	CPP01-1LT	Inquire	Inquire
Callus Proliferation Medium/ With macro- and micronutrients and vitamins.	Powder	10	L	CPP01-10LT	Inquire	Inquire
Embryo Production Medium/ With macro- and micronutrients and vitamins.	Powder	0.5	L	EPP01-0.5LT	Inquire	Inquire
Embryo Production Medium/ With macro- and micronutrients and vitamins.	Powder	1	L	EPP01-1LT	Inquire	Inquire
Embryo Production Medium/ With macro- and micronutrients and vitamins.	Powder	10	L	EPP01-10LT	Inquire	Inquire

DKW						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
DKW Medium/ With macro- and micronutrients	Powder	1	L	DKP01-1LT	Inquire	Inquire
DKW Medium/ With macro- and micronutrients	Powder	10	L	DKP01-10LT	Inquire	Inquire
DKW Medium/ With macro- and micronutrients	Powder	50	L	DKP01-50LT	Inquire	Inquire

E3P						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
E3P Medium with Vitamins/ Proprietary herbacious plant medium with macro- and micronutrients and vitamins	Powder	1	L	ETP01-1LT	Inquire	Inquire
E3P Medium with Vitamins/ Proprietary herbacious plant medium with macro- and micronutrients and vitamins	Powder	10	L	ETP01-10LT	Inquire	Inquire
E3P Medium with Vitamins/ Proprietary herbacious plant medium with macro- and micronutrients and vitamins	Powder	50	L	ETP01-50LT	Inquire	Inquire
E3P Medium with Vitamins and Adenine/ Proprietary herbacious plant medium with macro- and micronutrients, vitamins, and adenine hemisulfate	Powder	1	L	ETP02-1LT	Inquire	Inquire
E3P Medium with Vitamins and Adenine/ Proprietary herbacious plant medium with macro- and micronutrients, vitamins, and adenine hemisulfate	Powder	10	L	ETP02-10LT	Inquire	Inquire
E3P Medium with Vitamins and Adenine/ Proprietary herbacious plant medium with macro- and micronutrients, vitamins, and adenine hemisulfate	Powder	50	L	ETP02-50LT	Inquire	Inquire

Finer & Nagasawa						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Finer & Nagasawa Lite Modified Medium/ With MS micronutrients, Gamborg's vitamins, 1 g/L L-asparagine, and 10g/L sucrose	Powder	1	L	FNP01-1LT	Inquire	Inquire
Finer & Nagasawa Lite Modified Medium/ With MS micronutrients, Gamborg's vitamins, 1 g/L L-asparagine, and 10g/L sucrose	Powder	10	L	FNP01-10LT	Inquire	Inquire
Finer & Nagasawa Lite Modified Medium/ With MS micronutrients, Gamborg's vitamins, 1 g/L L-asparagine, and 10g/L sucrose	Powder	50	L	FNP01-50LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Gamborg's						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Gamborg's B-5 Long Medium/ With macro- and micronutrients (Gamborg et al., 1968) and enriched with 16 organic compounds	Powder	1	L	GBP06-1LT	Inquire	Inquire
Gamborg's B-5 Long Medium/ With macro- and micronutrients (Gamborg et al., 1968) and enriched with 16 organic compounds	Powder	10	L	GBP06-10LT	Inquire	Inquire
Gamborg's B-5 Long Medium/ With macro- and micronutrients (Gamborg et al., 1968) and enriched with 16 organic compounds	Powder	50	L	GBP06-50LT	Inquire	Inquire
Gamborg's B-5 Medium/ With macro- and micronutrients	Powder	1	L	GBP01-1LT	Inquire	Inquire
Gamborg's B-5 Medium/ With macro- and micronutrients	Powder	10	L	GBP01-10LT	Inquire	Inquire
Gamborg's B-5 Medium/ With macro- and micronutrients	Powder	50	L	GBP01-50LT	Inquire	Inquire
Gamborg's B-5 Medium with Sucrose/ With macro- and micronutrients and 20g/L sucrose	Powder	1	L	GBP07-1LT	Inquire	Inquire
Gamborg's B-5 Medium with Sucrose/ With macro- and micronutrients and 20g/L sucrose	Powder	10	L	GBP07-10LT	Inquire	Inquire
Gamborg's B-5 Medium with Sucrose/ With macro- and micronutrients and 20g/L sucrose	Powder	50	L	GBP07-50LT	Inquire	Inquire
Gamborg's B-5 Medium with Sucrose/ With macro- and micronutrients and 20g/L sucrose	Powder	50	L	GBP07-50LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Gamborg's B-5 Medium with Vitamins/ With macro- and micronutrients and vitamins.	Powder	1	L	GBP05-1LT	Inquire	Inquire
Gamborg's B-5 Medium with Vitamins/ With macro- and micronutrients and vitamins.	Powder	10	L	GBP05-10LT	Inquire	Inquire
Gamborg's B-5 Medium with Vitamins/ With macro- and micronutrients and vitamins.	Powder	50	L	GBP05-50LT	Inquire	Inquire
Gamborg's B-5 Medium with Vitamins and Sucrose/ With macro- and micronutrients, vitamins, and 20g/L sucrose.	Powder	1	L	GBP08-1LT	Inquire	Inquire
Gamborg's B-5 Medium with Vitamins and Sucrose/ With macro- and micronutrients, vitamins, and 20g/L sucrose.	Powder	10	L	GBP08-10LT	Inquire	Inquire
Gamborg's B-5 Medium with Vitamins and Sucrose/ With macro- and micronutrients, vitamins, and 20g/L sucrose.	Powder	50	L	GBP08-50LT	Inquire	Inquire
Gamborg's Vitamin Mixture/ 1000X Solution	Liquid	100	ML	GBL01-100ML	Inquire	Inquire

Heller's

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Heller's Medium/ With macro- and micronutrients	Powder	1	L	HEP01-1LT	Inquire	Inquire
Heller's Medium/ With macro- and micronutrients	Powder	10	L	HEP01-10LT	Inquire	Inquire
Heller's Medium/ With macro- and micronutrients	Powder	50	L	HEP01-50LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Hoagland's						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Hoagland's No. 2 Basal Salt Mixture/ With macro- and micronutrients	Powder	1	L	HOP01-1LT	Inquire	Inquire
Hoagland's No. 2 Basal Salt Mixture/ With macro- and micronutrients	Powder	10	L	HOP01-10LT	Inquire	Inquire
Hoagland's No. 2 Basal Salt Mixture/ With macro- and micronutrients	Powder	50	L	HOP01-50LT	Inquire	Inquire
Hoagland's No. 2 Basal Salt Mixture without Nitrogen/ With macro- and micronutrients; Without nitrogen	Powder	1	L	HOP03-1LT	Inquire	Inquire
Hoagland's No. 2 Basal Salt Mixture without Nitrogen/ With macro- and micronutrients; Without nitrogen	Powder	10	L	HOP03-10LT	Inquire	Inquire
Hoagland's No. 2 Basal Salt Mixture without Nitrogen/ With macro- and micronutrients; Without nitrogen	Powder	50	L	HOP03-50LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Kao & Michayluk						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Kao & Michayluk/ With macro- and micronutrients	Powder	1	L	KMP01-1LT	Inquire	Inquire
Kao & Michayluk/ With macro- and micronutrients	Powder	10	L	KMP01-10LT	Inquire	Inquire
Kao & Michayluk/ With macro- and micronutrients	Powder	50	L	KMP01-50LT	Inquire	Inquire
Kao & Michayluk Vitamin Mixture/ 100X Solution	Liquid	100	ML	KML01-100ML	Inquire	Inquire
Kao & Michayluk with Vitamins/ With macro- and micronutrients and vitamins.	Powder	1	L	KMP02-1LT	Inquire	Inquire
Kao & Michayluk with Vitamins/ With macro- and micronutrients and vitamins.	Powder	10	L	KMP02-10LT	Inquire	Inquire
Kao & Michayluk with Vitamins/ With macro- and micronutrients and vitamins.	Powder	50	L	KMP02-50LT	Inquire	Inquire

Lindemann						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Lindemann Orchid Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	1	L	LOP01-1LT	Inquire	Inquire
Lindemann Orchid Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	1	L	LOP01-1LT	Inquire	Inquire
Lindemann Orchid Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	10	L	LOP01-10LT	Inquire	Inquire
Lindemann Orchid Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	50	L	LOP01-50LT	Inquire	Inquire

Linsmaier & Skoog						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Linsmaier & Skoog/ With macro- and micronutrients and vitamins	Powder	1	L	LSP01-1LT	Inquire	Inquire
Linsmaier & Skoog/ With macro- and micronutrients and vitamins	Powder	10	L	LSP01-10LT	Inquire	Inquire
Linsmaier & Skoog/ With macro- and micronutrients and vitamins	Powder	50	L	LSP01-50LT	Inquire	Inquire
Linsmaier & Skoog with Buffer/ With macro- and micronutrients; pH buffered to 5.7 +/-0.2 and vitamins	Powder	1	L	LSP03-1LT	Inquire	Inquire
Linsmaier & Skoog with Buffer/ With macro- and micronutrients; pH buffered to 5.7 +/-0.2 and vitamins	Powder	10	L	LSP03-10LT	Inquire	Inquire
Linsmaier & Skoog with Buffer/ With macro- and micronutrients; pH buffered to 5.7 +/-0.2 and vitamins	Powder	50	L	LSP03-50LT	Inquire	Inquire
Linsmaier & Skoog with Buffer and Sucrose/ With macro- and micronutrients and 30 g/L sucrose; pH buffered to 5.7 +/-0.2 and vitamins	Powder	1	L	LSP04-1LT	Inquire	Inquire
Linsmaier & Skoog with Buffer and Sucrose/ With macro- and micronutrients and 30 g/L sucrose; pH buffered to 5.7 +/-0.2 and vitamins	Powder	1	L	LSP04-1LT	Inquire	Inquire

Murashige & Skoog						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Murashige & Skoog Basal Salts/ With macro- and micronutrients	Powder	1	L	MSP01-1LT	Inquire	Inquire
Murashige & Skoog Basal Salts/ With macro- and micronutrients	Powder	10	L	MSP01-10LT	Inquire	Inquire
Murashige & Skoog Basal Salts/ With macro- and micronutrients	Powder	50	L	MSP01-50LT	Inquire	Inquire
Murashige & Skoog Basal Salts with Sucrose/ With macro- and micronutrients and 30g/l sucrose	Powder	1	L	MSP13-1LT	Inquire	Inquire
Murashige & Skoog Basal Salts with Sucrose/ With macro- and micronutrients and 30g/l sucrose	Powder	10	L	MSP13-10LT	Inquire	Inquire
Murashige & Skoog Basal Salts with Sucrose/ With macro- and micronutrients and 30g/l sucrose	Powder	50	L	MSP13-50LT	Inquire	Inquire
Murashige & Skoog Macronutrients/ Macronutrients only	Powder	1	L	MSP17-1LT	Inquire	Inquire
Murashige & Skoog Macronutrients/ Macronutrients only	Powder	10	L	MSP17-10LT	Inquire	Inquire
Murashige & Skoog Macronutrients/ Macronutrients only	Powder	50	L	MSP17-50LT	Inquire	Inquire
Murashige & Skoog Micronutrient Stock/ Miconutrients only; Each liter of powder makes one liter of a 100X stock solution	Powder	1	L	MSP18-1LT	Inquire	Inquire
Murashige & Skoog Micronutrient Stock/ Miconutrients only; Each liter of powder makes one liter of a 100X stock solution	Powder	10	L	MSP18-10LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

<p>Murashige & Skoog Micronutrient Stock/ Miconutrients only; Each liter of powder makes one liter of a 100X stock solution</p>	Powder	50	L	MSP18-50LT	Inquire	Inquire
<p>Murashige & Skoog Modified with Adenine, Inositol, and Thiamine/ With macro- and micronutrients, 80mg/L adenine sulfate, 100mg/L myo-inositol, and 0.40mg/L thiamine HCl</p>	Powder	1	L	MSP03-1LT	Inquire	Inquire
<p>Murashige & Skoog Modified with Adenine, Inositol, and Thiamine/ With macro- and micronutrients, 80mg/L adenine sulfate, 100mg/L myo-inositol, and 0.40mg/L thiamine HCl</p>	Powder	10	L	MSP03-10LT	Inquire	Inquire
<p>Murashige & Skoog Modified with Adenine, Inositol, and Thiamine/ With macro- and micronutrients, 80mg/L adenine sulfate, 100mg/L myo-inositol, and 0.40mg/L thiamine HCl</p>	Powder	50	L	MSP03-50LT	Inquire	Inquire
<p>Murashige & Skoog Multiplication Medium C (MSMC)/ With macro- and micronutrients, vitamins, sucrose, sodium phosphate, NAA, adenine sulfate, and kinetin</p>	Powder	1	L	MSP12-1LT	Inquire	Inquire
<p>Murashige & Skoog Multiplication Medium C (MSMC)/ With macro- and micronutrients, vitamins, sucrose, sodium phosphate, NAA, adenine sulfate, and kinetin</p>	Powder	10	L	MSP12-10LT	Inquire	Inquire
<p>Murashige & Skoog Multiplication Medium C (MSMC)/ With macro- and micronutrients, vitamins, sucrose, sodium phosphate, NAA, adenine sulfate, and kinetin</p>	Powder	50	L	MSP12-50LT	Inquire	Inquire
<p>Murashige & Skoog Vitamin Mixture/ 1000X Solution</p>	Liquid	100	ML	MVL01-100ML	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Murashige & Skoog with Gamborg's Vitamins/ With macro- and micronutrients and Gamborg's vitamins	Powder	1	L	MSP06-1LT	Inquire	Inquire
Murashige & Skoog with Gamborg's Vitamins/ With macro- and micronutrients and Gamborg's vitamins	Powder	10	L	MSP06-10LT	Inquire	Inquire
Murashige & Skoog with Gamborg's Vitamins/ With macro- and micronutrients and Gamborg's vitamins	Powder	50	L	MSP06-50LT	Inquire	Inquire
Murashige & Skoog with Reduced Macronutrients and Iron/ With macro- and micronutrients; modified with macronutrients and iron concentrations cut in half	Powder	1	L	MSP08-1LT	Inquire	Inquire
Murashige & Skoog with Reduced Macronutrients and Iron/ With macro- and micronutrients; modified with macronutrients and iron concentrations cut in half	Powder	10	L	MSP08-10LT	Inquire	Inquire
Murashige & Skoog with Reduced Macronutrients and Iron/ With macro- and micronutrients; modified with macronutrients and iron concentrations cut in half	Powder	50	L	MSP08-50LT	Inquire	Inquire
Murashige & Skoog with Reduced NH4NO3 and KNO3/ With macro- and micronutrients; modified with NH4NO3 and KNO3 concentrations cut in half	Powder	1	L	MSP16-1LT	Inquire	Inquire
Murashige & Skoog with Reduced NH4NO3 and KNO3/ With macro- and micronutrients; modified with NH4NO3 and KNO3 concentrations cut in half	Powder	10	L	MSP16-10LT	Inquire	Inquire
Murashige & Skoog with Reduced NH4NO3 and KNO3/ With macro- and micronutrients; modified with NH4NO3 and KNO3 concentrations cut in half	Powder	50	L	MSP16-50LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

<p>Murashige & Skoog with Reduced NH₄NO₃, KNO₃, and CaCl₂/</p> <p>With macro- and micronutrients; modified with NH₄NO₃, KNO₃, and CaCl₂ concentrations cut in half</p>	Powder	1	L	MSP04-1LT	Inquire	Inquire
<p>Murashige & Skoog with Reduced NH₄NO₃, KNO₃, and CaCl₂/</p> <p>With macro- and micronutrients; modified with NH₄NO₃, KNO₃, and CaCl₂ concentrations cut in half</p>	Powder	10	L	MSP04-10LT	Inquire	Inquire
<p>Murashige & Skoog with Reduced NH₄NO₃, KNO₃, and CaCl₂/</p> <p>With macro- and micronutrients; modified with NH₄NO₃, KNO₃, and CaCl₂ concentrations cut in half</p>	Powder	50	L	MSP04-50LT	Inquire	Inquire
<p>Murashige & Skoog with Vitamins/</p> <p>With macro- and micronutrients and vitamins</p>	Powder	1	L	MSP09-1LT	Inquire	Inquire
<p>Murashige & Skoog with Vitamins/</p> <p>With macro- and micronutrients and vitamins</p>	Powder	10	L	MSP09-10LT	Inquire	Inquire
<p>Murashige & Skoog with Vitamins/</p> <p>With macro- and micronutrients and vitamins</p>	Powder	50	L	MSP09-50LT	Inquire	Inquire
<p>Murashige & Skoog with Vitamins and Sucrose/</p> <p>With macro- and micronutrients, vitamins, and 30g/L sucrose</p>	Powder	1	L	MSP14-1LT	Inquire	Inquire
<p>Murashige & Skoog with Vitamins and Sucrose/</p> <p>With macro- and micronutrients, vitamins, and 30g/L sucrose</p>	Powder	10	L	MSP14-10LT	Inquire	Inquire
<p>Murashige & Skoog with Vitamins and Sucrose/</p> <p>With macro- and micronutrients, vitamins, and 30g/L sucrose</p>	Powder	50	L	MSP14-50LT	Inquire	Inquire
<p>Murashige & Skoog with Vitamins; without Glycine/</p> <p>With macro- and micronutrients and vitamins; without glycine</p>	Powder	1	L	MSP02-1LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Murashige & Skoog with Vitamins; without Glycine/ With macro- and micronutrients and vitamins; without glycine	Powder	10	L	MSP02-10LT	Inquire	Inquire
Murashige & Skoog with Vitamins; without Glycine/ With macro- and micronutrients and vitamins; without glycine	Powder	50	L	MSP02-50LT	Inquire	Inquire
Murashige & Skoog without NH₄NO₃/ With macro- and micronutrients; without NH ₄ NO ₃	Powder	1	L	MSP05-1LT	Inquire	Inquire
Murashige & Skoog without NH₄NO₃/ With macro- and micronutrients; without NH ₄ NO ₃	Powder	10	L	MSP05-10LT	Inquire	Inquire
Murashige & Skoog without NH₄NO₃/ With macro- and micronutrients; without NH ₄ NO ₃	Powder	50	L	MSP05-50LT	Inquire	Inquire
Murashige & Skoog without Nitrates/ With nitrogen-free macro- and micronutrients and vitamins.	Powder	10	L	MSP07-10LT	Inquire	Inquire
Murashige & Skoog without Nitrates/ With nitrogen-free macro- and micronutrients and vitamins.	Powder	50	L	MSP07-50LT	Inquire	Inquire
Murashige & Skoog without Nitrogen/ With nitrogen-free macro- and micronutrients. Without vitamins.	Powder	1	L	MSP21-1LT	Inquire	Inquire
Murashige & Skoog without Nitrogen/ With nitrogen-free macro- and micronutrients. Without vitamins.	Powder	10	L	MSP21-10LT	Inquire	Inquire
Murashige & Skoog without Nitrogen/ With nitrogen-free macro- and micronutrients. Without vitamins.	Powder	50	L	MSP21-50LT	Inquire	Inquire
Murashige & Skoog without Phosphate/ With macro- and micronutrients; without phosphate	Powder	1	L	MSP11-1LT	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Murashige & Skoog without Phosphate/ With macro- and micronutrients; without phosphate	Powder	10	L	MSP11-10LT	Inquire	Inquire
Murashige & Skoog without Phosphate/ With macro- and micronutrients; without phosphate	Powder	50	L	MSP11-50LT	Inquire	Inquire

Nitsch & Nitsch

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Nitsch & Nitsch/ With macro- and micronutrients	Powder	1	L	NNP01-1LT	Inquire	Inquire
Nitsch & Nitsch/ With macro- and micronutrients	Powder	10	L	NNP01-10LT	Inquire	Inquire
Nitsch & Nitsch/ With macro- and micronutrients	Powder	50	L	NNP01-50LT	Inquire	Inquire
Nitsch & Nitsch with Vitamins/ With macro- and micronutrients and vitamins	Powder	1	L	NNP03-1LT	Inquire	Inquire
Nitsch & Nitsch with Vitamins/ With macro- and micronutrients and vitamins	Powder	10	L	NNP03-10LT	Inquire	Inquire
Nitsch & Nitsch with Vitamins/ With macro- and micronutrients and vitamins	Powder	50	L	NNP03-50LT	Inquire	Inquire

Orchid Maintenance

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Orchid Maintenance Medium/ Without agar and charcoal	Powder	1	L	OMP01-1LT	Inquire	Inquire
Orchid Maintenance Medium/ Without agar and charcoal	Powder	10	L	OMP01-10LT	Inquire	Inquire
Orchid Maintenance Medium/ Without agar and charcoal	Powder	50	L	OMP01-50LT	Inquire	Inquire

Quoirin & Lepoire						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Quoirin & Lepoire Basal Salt Mixture/ With macro- and micronutrients	Powder	1	L	QLP01-1LT	Inquire	Inquire
Quoirin & Lepoire Basal Salt Mixture/ With macro- and micronutrients	Powder	10	L	QLP01-10LT	Inquire	Inquire
Quoirin & Lepoire Basal Salt Mixture/ With macro- and micronutrients	Powder	50	L	QLP01-50LT	Inquire	Inquire

Shenk & Hildebrandt						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Shenk & Hildebrandt Basal Salts/ With macro- and micronutrients (Shenk & Hildebrandt , 1972)	Powder	1	L	SHP01-1LT	Inquire	Inquire
Shenk & Hildebrandt Basal Salts/ With macro- and micronutrients (Shenk & Hildebrandt , 1972)	Powder	10	L	SHP01-10LT	Inquire	Inquire
Shenk & Hildebrandt Basal Salts/ With macro- and micronutrients (Shenk & Hildebrandt , 1972)	Powder	50	L	SHP01-50LT	Inquire	Inquire

Vacin & Went						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Vacin & Went Orchid Medium/ With macro- and micronutrients (Vacin & Went, 1949)	Powder	1	L	VWP01-1LT	Inquire	Inquire
Vacin & Went Orchid Medium/ With macro- and micronutrients (Vacin & Went, 1949)	Powder	10	L	VWP01-10LT	Inquire	Inquire
Vacin & Went Orchid Medium/ With macro- and micronutrients (Vacin & Went, 1949)	Powder	50	L	VWP01-50LT	Inquire	Inquire

Vitamin Mixtures						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Gamborg's Vitamin Mixture/ 1000X Solution	Liquid	100	ML	GBL01-100ML	Inquire	Inquire
Kao & Michayluk Vitamin Mixture/ 100X Solution	Liquid	100	ML	KML01-100ML	Inquire	Inquire
Murashige & Skoog Vitamin Mixture/ 1000X Solution	Liquid	100	ML	MVL01-100ML	Inquire	Inquire

Product Catalog

Plant Tissue Culture Tissue Culture Media

Westvaco						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Westvaco WV3 Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	1	L	WVP01-1LT	Inquire	Inquire
Westvaco WV3 Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	10	L	WVP01-10LT	Inquire	Inquire
Westvaco WV3 Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	50	L	WVP01-50LT	Inquire	Inquire
Westvaco WV5 Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	1	L	WVP02-1LT	Inquire	Inquire
Westvaco WV5 Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	10	L	WVP02-10LT	Inquire	Inquire
Westvaco WV5 Medium with Vitamins/ With macro- and micronutrients and vitamins	Powder	50	L	WVP02-50LT	Inquire	Inquire

Woody Plant						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Lloyd & McKown's Woody Plant Medium/ With macro- and micronutrients	Powder	1	L	WPP01-1LT	Inquire	Inquire
Lloyd & McKown's Woody Plant Medium/ With macro- and micronutrients	Powder	10	L	WPP01-10LT	Inquire	Inquire
Lloyd & McKown's Woody Plant Medium/ With macro- and micronutrients	Powder	50	L	WPP01-50LT	Inquire	Inquire

Plant Tissue Culture
Vessels and Supplies

Product Catalog

Plant Tissue Culture Vessels and Supplies

Caisson Brand						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
MK-5 (Case of 100)/ 3"x3"x4" Polycarbonate Vessel with Cover; Case of 100; Compare to Magenta GA-7		100	PCS	MK5-100PC	Inquire	Inquire
MK-5 (Case of 20)/ 3"x3"x4" Polycarbonate Vessel with Cover; Case of 20; Compare to Magenta GA-7		20	PCS	MK5-20PC	Inquire	Inquire
Vessel Wrap (Sealing Tape)/ 80 gauge; 1.5" width x 1000' length. 1 Roll		1	PC	VWS01-1PC	Inquire	Inquire
Vessel Wrap (Sealing Tape)/ 80 gauge; 1.5" width x 1000' length per roll. 10 Rolls		10	PC	VWS01-10PC	Inquire	Inquire

Liquid Lab						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Liquid Lab Vessel/ 10.5"x 4"x 4"; Polycarbonate; Comes with microporous patch and a septaport		1	PCS	SSV01-1PC	Inquire	Inquire
Liquid Lab Vessel/ 10.5"x 4"x 4"; Polycarbonate; Comes with microporous patch and a septaport		1	PCS	SSV01-1PC	Inquire	Inquire

Product Catalog

Plant Tissue Culture Vessels and Supplies

Magenta						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
B-Caps/ Caps for baby food jars; Polypropylene; Pack of 100		100	PCS	BCAP-100PC	Inquire	Inquire
Coupler/ Coupler for GA-7, GA-7-3, and MK-5; Allows another vessel to be used in place of a lid thus doubling height of growing space.		100	PCS	GA7C-100PC	Inquire	Inquire
GA-7 (Case of 100)/ 3"x3"x4" Polycarbonate Vessel with Cover; Case of 100		100	PCS	GA7-100PC	Inquire	Inquire
GA-7-3 (Case of 100)/ 3"x3"x3" Polycarbonate Vessel with Cover; Case of 100		100	PCS	GA73-100PC	Inquire	Inquire
Test Tube Cap/ Autoclavable caps for 25mm test tubes		100	PCS	MTTC-100PC	Inquire	Inquire

Petri Dishes						
NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
100 x 15mm Petri Dishes/ 100 x 15mm; Sterile; Disposable; Case of 500		500	PCS	PET01-500PC	Inquire	Inquire