



Product Catalog 2017

Antibiotics
Anti-Bacterial

Product Catalog Antibiotics Anti-Bacterial

Ampicillin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Ampicillin, Sodium Salt/	Powder	5g	Bottle	A043-5GM	Inquire	Inquire
Ampicillin, Sodium Salt/	Powder	25g	Bottle	A043-25GM	Inquire	Inquire
Ampicillin, Sodium Salt/	Powder	100g	Bottle	A043-100GM	Inquire	Inquire

Carbenicillin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Carbenicillin, Disodium Salt/ USP Grade; CAS Number: 4800-94-6; Formula: C ₁₇ H ₁₆ N ₂ O ₆ SNa ₂ ; Formula Weight: 422.36. Store at 2° - 8°C.	Powder	5g	Bottle	C033-5GM	Inquire	Inquire
Carbenicillin, Disodium Salt/ USP Grade; CAS Number: 4800-94-6; Formula: C ₁₇ H ₁₆ N ₂ O ₆ SNa ₂ ; Formula Weight: 422.36. Store at 2° - 8°C.	Powder	10g	Bottle	C033-10GM	Inquire	Inquire
Carbenicillin, Disodium Salt/ USP Grade; CAS Number: 4800-94-6; Formula: C ₁₇ H ₁₆ N ₂ O ₆ SNa ₂ ; Formula Weight: 422.36. Store at 2° - 8°C.	Powder	25g	Bottle	C033-25GM	Inquire	Inquire
Carbenicillin, Disodium Salt/ USP Grade; CAS Number: 4800-94-6; Formula: C ₁₇ H ₁₆ N ₂ O ₆ SNa ₂ ; Formula Weight: 422.36. Store at 2° - 8°C.	Powder	100g	Bottle	C033-100GM	Inquire	Inquire

Product Catalog Antibiotics Anti-Bacterial

Cefotaxime

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Cefotaxime, Sodium Salt/ USP Grade; CAS Number: 64485-93-4; Formula: C ₁₆ H ₁₆ N ₅ NaO ₇ S ₂ ; Formula Weight: 477.45. Store at 2° - 8°C.	Powder	1g	Bottle	C032-1GM	Inquire	Inquire
Cefotaxime, Sodium Salt/ USP Grade; CAS Number: 64485-93-4; Formula: C ₁₆ H ₁₆ N ₅ NaO ₇ S ₂ ; Formula Weight: 477.45. Store at 2° - 8°C.	Powder	5g	Bottle	C032-5GM	Inquire	Inquire
Cefotaxime, Sodium Salt/ USP Grade; CAS Number: 64485-93-4; Formula: C ₁₆ H ₁₆ N ₅ NaO ₇ S ₂ ; Formula Weight: 477.45. Store at 2° - 8°C.	Powder	10g	Bottle	C032-10GM	Inquire	Inquire
Cefotaxime, Sodium Salt/ USP Grade; CAS Number: 64485-93-4; Formula: C ₁₆ H ₁₆ N ₅ NaO ₇ S ₂ ; Formula Weight: 477.45. Store at 2° - 8°C.	Powder	25g	Bottle	C032-25GM	Inquire	Inquire
Cefotaxime, Sodium Salt/ USP Grade; CAS Number: 64485-93-4; Formula: C ₁₆ H ₁₆ N ₅ NaO ₇ S ₂ ; Formula Weight: 477.45. Store at 2° - 8°C.	Powder	100g	Bottle	C032-100GM	Inquire	Inquire

Product Catalog Antibiotics Anti-Bacterial

Gentamicin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Gentamicin Solution (50mg/ml)/	Liquid	10ml	Bottle	ABL03-10ML	Inquire	Inquire
Gentamicin Solution (50mg/ml)/	Liquid	6 x 10ml	Bottles	ABL03-6X10ML	Inquire	Inquire
Gentamicin Solution (10mg/ml)/	Liquid	10ml	Bottle	ABL05-10ML	Inquire	Inquire
Gentamicin Solution (10mg/ml)/	Liquid	6 x 10ml	Bottles	ABL05-6X10ML	Inquire	Inquire
Gentamicin Sulfate/	Powder	1g	Bottle	G020-1GM	Inquire	Inquire
Gentamicin Sulfate/	Powder	5g	Bottle	G020-5GM	Inquire	Inquire
Gentamicin Sulfate/	Powder	10g	Bottle	G020-10GM	Inquire	Inquire
Gentamicin Sulfate/	Powder	25g	Bottle	G020-25GM	Inquire	Inquire
Gentamicin Sulfate/	Powder	100g	Bottle	G020-100GM	Inquire	Inquire

Penicillin G

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Penicillin G, Sodium Salt/	Powder	25g	Bottle	P056-25GM	Inquire	Inquire
Penicillin G, Sodium Salt/	Powder	100g	Bottle	P056-100GM	Inquire	Inquire

Penicillin-Streptomycin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
<p>Penicillin-Streptomycin Solution 100X/ This sterile-filtered solution contains 10,000 units/ml of penicillin and 10,000 µg/ml streptomycin in 0.9% saline. Penicillin-Streptomycin solution is effective against gram-positive and gram-negative bacteria. Suitable for use in cell culture. Store at -5° to -20°C.</p>	Liquid	100ml	Bottle	PSL01-100ML	Inquire	Inquire
<p>Penicillin-Streptomycin Solution 100X/ This sterile-filtered solution contains 10,000 units/ml of penicillin and 10,000 µg/ml streptomycin in 0.9% saline. Penicillin-Streptomycin solution is effective against gram-positive and gram-negative bacteria. Suitable for use in cell culture. Store at -5° to -20°C.</p>	Liquid	6 x 100ml	Bottles	PSL01-6X100ML	Inquire	Inquire
<p>Penicillin-Streptomycin-Glutamine Solution 100X/ This sterile-filtered solution contains 10,000 units/ml of penicillin, 10,000 µg/ml streptomycin and 29.2 mg/ml L-glutamine in 0.9% saline. Penicillin-Streptomycin-Glutamine solution is effective against gram-positive and gram-negative bacteria. Suitable for use in cell culture. Store at -5° to -20°C.</p>	Liquid	100ml	Bottle	PSL02-100ML	Inquire	Inquire

Product Catalog Antibiotics Anti-Bacterial

<p>Penicillin-Streptomycin-Glutamine Solution 100X/</p> <p>This sterile-filtered solution contains 10,000 units/ml of penicillin, 10,000 µg/ml streptomycin and 29.2 mg/ml L-glutamine in 0.9% saline. Penicillin-Streptomycin-Glutamine solution is effective against gram-positive and gram-negative bacteria. Suitable for use in cell culture. Store at -5° to -20°C.</p>	Liquid	6 x 100ml	Bottles	PSL02-6X100ML	Inquire	Inquire
--	--------	--------------	---------	---------------	---------	---------

Puromycin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Puromycin Dihydrochloride/	Powder	100mg	Bottle	P067-100MG	Inquire	Inquire
Puromycin Dihydrochloride/	Powder	500mg	Bottle	P067-500MG	Inquire	Inquire
Puromycin Dihydrochloride/	Powder	1g	Bottle	P067-1GM	Inquire	Inquire

Streptomycin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Streptomycin Sulfate/	Powder	25g	Bottle	S041-25GM	Inquire	Inquire
Streptomycin Sulfate/	Powder	100g	Bottle	S041-100GM	Inquire	Inquire

Product Catalog Antibiotics Anti-Bacterial

Timentin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Timentin/ 50X (50mg/ml) Timentin (15:1 Ticarcillin Disodium Salt : Potassium Clavulanate Mixture) solution. Sterile filtered. Store at -20 to -5°C.	Liquid	10ml	Bottle	ABL07-10ML	Inquire	Inquire
Timentin/ 50X (50mg/ml) Timentin (15:1 Ticarcillin Disodium Salt : Potassium Clavulanate Mixture) solution. Sterile filtered. Store at -20 to -5°C.	Liquid	50ml	Bottle	ABL07-50ML	Inquire	Inquire
Timentin/ 50X (50mg/ml) Timentin (15:1 Ticarcillin Disodium Salt : Potassium Clavulanate Mixture) solution. Sterile filtered. Store at -20 to -5°C.	Liquid	100ml	Bottle	ABL07-100ML	Inquire	Inquire
Timentin/ 100X (100mg/ml) Timentin (15:1 Ticarcillin Disodium Salt : Potassium Clavulanate Mixture) solution. Sterile filtered. Store at -20 to -5°C.	Liquid	10ml	Bottle	ABL08-10ML	Inquire	Inquire
Timentin/ 100X (100mg/ml) Timentin (15:1 Ticarcillin Disodium Salt : Potassium Clavulanate Mixture) solution. Sterile filtered. Store at -20 to -5°C.	Liquid	50ml	Bottle	ABL08-50ML	Inquire	Inquire
Timentin/ 100X (100mg/ml) Timentin (15:1 Ticarcillin Disodium Salt : Potassium Clavulanate Mixture) solution. Sterile filtered. Store at -20 to -5°C.	liquid	100ml	Bottle	ABL08-100ML	Inquire	Inquire
Timentin/ 15:1 Ticarcillin Disodium Salt / Potassium Clavulanate Mixture; CAS Numbers: 4697-14-7 and 61177-45-5. Store at 2° to 8°C.	Powder	1g	Bottle	T034-1GM	Inquire	Inquire
Timentin/ 15:1 Ticarcillin Disodium Salt / Potassium Clavulanate Mixture; CAS Numbers: 4697-14-7 and 61177-45-5. Store at 2° to 8°C.	Powder	5g	Bottle	T034-5GM	Inquire	Inquire

Product Catalog Antibiotics Anti-Bacterial

<p>Timentin/ 15:1 Ticarcillin Disodium Salt / Potassium Clavulanate Mixture; CAS Numbers: 4697-14-7 and 61177-45-5. Store at 2° to 8°C.</p>	Powder	25g	Bottle	T034-25GM	Inquire	Inquire
<p>Timentin/ 15:1 Ticarcillin Disodium Salt / Potassium Clavulanate Mixture; CAS Numbers: 4697-14-7 and 61177-45-5. Store at 2° to 8°C.</p>	Powder	100g	Bottle	T034-100GM	Inquire	Inquire

Vancomycin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
<p>Vancomycin HCl/ Used in combination with Cefotaxime to kill agrobacterium. CAS Number: 1404-93-9; Formula: C₆₆H₇₅Cl₂N₉O₂₄ . HCl; Formula Weight: 1485.71; Assay (dried basis): >900 ug/mg. Store below -20°C.</p>	Powder	1g	Bottle	V007-1GM	Inquire	Inquire
<p>Vancomycin HCl/ Used in combination with Cefotaxime to kill agrobacterium. CAS Number: 1404-93-9; Formula: C₆₆H₇₅Cl₂N₉O₂₄ . HCl; Formula Weight: 1485.71; Assay (dried basis): >900 ug/mg. Store below -20°C.</p>	Powder	5g	Bottle	V007-5GM	Inquire	Inquire
<p>Vancomycin HCl/ Used in combination with Cefotaxime to kill agrobacterium. CAS Number: 1404-93-9; Formula: C₆₆H₇₅Cl₂N₉O₂₄ . HCl; Formula Weight: 1485.71; Assay (dried basis): >900 ug/mg. Store below -20°C.</p>	Powder	10g	Bottle	V007-10GM	Inquire	Inquire
<p>Vancomycin HCl/ Used in combination with Cefotaxime to kill agrobacterium. CAS Number: 1404-93-9; Formula: C₆₆H₇₅Cl₂N₉O₂₄ . HCl; Formula Weight: 1485.71; Assay (dried basis): >900 ug/mg. Store below -20°C.</p>	Powder	25g	Bottle	V007-25GM	Inquire	Inquire

Antibiotics
Anti-Mycotic

Product Catalog

Antibiotics Anti-Mycotic

Amphotericin B

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Amphotericin B/	Powder	500mg	Bottle	A039-500MG	Inquire	Inquire
Amphotericin B/	Powder	1g	Bottle	A039-1GM	Inquire	Inquire
Amphotericin B/	Powder	5g	Bottle	A039-5GM	Inquire	Inquire
Amphotericin B Solution/	Liquid	100ml	Bottle	ABL01-100ML	Inquire	Inquire
Amphotericin B Solution/	liquid	6 x 100ml	Bottles	ABL01-6X100ML	Inquire	Inquire

Nystatin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Nystatin Suspension/	Liquid	100	ml	ABL04-100ML	Inquire	Inquire
Nystatin Suspension/	Liquid	6 x100ml	Bottles	ABL04-6X100ML	Inquire	Inquire
Nystatin/	Powder	5g	Bottle	N011-5GM	Inquire	Inquire
Nystatin/	Powder	10g	Bottle	N011-10GM	Inquire	Inquire
Nystatin/	Powder	25g	Bottle	N011-25GM	Inquire	Inquire

Antibiotics
Anti-Viral

Rifampicin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Rifampicin/ CAS Number: 13292-46-1; Formula: C ₄₃ H ₅₈ N ₄ O ₁₂ ; Formula Weight: 822. 94; Store below 0°C.	Powder	1g	Bottle	R008-1GM	Inquire	Inquire
Rifampicin/ CAS Number: 13292-46-1; Formula: C ₄₃ H ₅₈ N ₄ O ₁₂ ; Formula Weight: 822. 94; Store below 0°C.	Powder	5g	Bottle	R008-5GM	Inquire	Inquire
Rifampicin/ CAS Number: 13292-46-1; Formula: C ₄₃ H ₅₈ N ₄ O ₁₂ ; Formula Weight: 822. 94; Store below 0°C.	Powder	10g	Bottle	R008-10GM	Inquire	Inquire
Rifampicin/ CAS Number: 13292-46-1; Formula: C ₄₃ H ₅₈ N ₄ O ₁₂ ; Formula Weight: 822. 94; Store below 0°C.	Powder	25g	Bottle	R008-25GM	Inquire	Inquire

Antibiotics Broad Range

Product Catalog Antibiotics Broad Range

Antibiotic-Antimycotic

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Antibiotic-Antimycotic Solution (100X)/	Liquid	100ml	Bottle	ABL02-100ML	Inquire	Inquire
Antibiotic-Antimycotic Solution (100X)/	Liquid	6 x 100ml	Bottles	ABL02-6X100ML	Inquire	Inquire

Plant Preservative Mixture (PPM)

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Plant Preservative Mixture (PPM™)/	Liquid	10ml	Bottle	PPL33-10ML	Inquire	Inquire
Plant Preservative Mixture (PPM™)/	Liquid	30ml	Bottle	PPL33-30ML	Inquire	Inquire
Plant Preservative Mixture (PPM™)/	Liquid	100ml	Bottle	PPL33-100ML	Inquire	Inquire
Plant Preservative Mixture (PPM™)/	Liquid	500ml	Bottle	PPL33-500ML	Inquire	Inquire
Plant Preservative Mixture (PPM™)/	Liquid	1000ml	Bottle	PPL33-1000ML	Inquire	Inquire

Antibiotics Selection

Bialaphos

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
<p>Bialaphos Sodium Salt/ Bialaphos is a tripeptide antibiotic with herbicidal properties used to select the bar gene. Bialaphos works great on monocots such as wheat and corn. CAS Number 7108-99-2; Formula: C₁₁H₂₁N₃NaO₆PNa; Formula Weight: 345.26. Store at -20 °C.</p>	Powder	100mg	Bottle	B009-100MG	Inquire	Inquire
<p>Bialaphos Sodium Salt/ Bialaphos is a tripeptide antibiotic with herbicidal properties used to select the bar gene. Bialaphos works great on monocots such as wheat and corn. CAS Number 7108-99-2; Formula: C₁₁H₂₁N₃NaO₆PNa; Formula Weight: 345.26. Store at -20 °C.</p>	Powder	250mg	Bottle	B009-250MG	Inquire	Inquire
<p>Bialaphos Sodium Salt/ Bialaphos is a tripeptide antibiotic with herbicidal properties used to select the bar gene. Bialaphos works great on monocots such as wheat and corn. CAS Number 7108-99-2; Formula: C₁₁H₂₁N₃NaO₆PNa; Formula Weight: 345.26. Store at -20 °C.</p>	Powder	500mg	Bottle	B009-500MG	Inquire	Inquire
<p>Bialaphos Sodium Salt/ Bialaphos is a tripeptide antibiotic with herbicidal properties used to select the bar gene. Bialaphos works great on monocots such as wheat and corn. CAS Number 7108-99-2; Formula: C₁₁H₂₁N₃NaO₆PNa; Formula Weight: 345.26. Store at -20 °C.</p>	Powder	1g	Bottle	B009-1GM	Inquire	Inquire

Product Catalog Antibiotics Selection

G-418

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
G418 Sulfate Solution/ G418 Sulfate Solution, 50 mg/ml. Store at 2° to 8°C. Ships with ice packs.	Liquid	20ml	Bottle	ABL06-20ML	Inquire	Inquire
G418 Sulfate Solution/ G418 Sulfate Solution, 50 mg/ml. Store at 2° to 8°C. Ships with ice packs.	Liquid	50ml	Bottle	ABL06-50ML	Inquire	Inquire
G418 Sulfate Solution/ G418 Sulfate Solution, 50 mg/ml. Store at 2° to 8°C. Ships with ice packs.	Liquid	100ml	Bottle	ABL06-100ML	Inquire	Inquire
G418 disulfate/	Powder	250mg	Bottle	G030-250MG	Inquire	Inquire
G418 disulfate/	Powder	1g	Bottle	G030-1GM	Inquire	Inquire
G418 disulfate/	Powder	5g	Bottle	G030-5GM	Inquire	Inquire
G418 disulfate/	Powder	10g	Bottle	G030-10GM	Inquire	Inquire

Gentamicin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Gentamicin Solution (50mg/ml)/	Liquid	10ml	Bottle	ABL03-10ML	Inquire	Inquire
Gentamicin Solution (50mg/ml)/	Liquid	6 x 10ml	Bottles	ABL03-6X10ML	Inquire	Inquire
Gentamicin Sulfate/	Powder	1g	Bottle	G020-1GM	Inquire	Inquire
Gentamicin Sulfate/	Powder	5g	Bottle	G020-5GM	Inquire	Inquire
Gentamicin Sulfate/	Powder	10g	Bottle	G020-10GM	Inquire	Inquire
Gentamicin Sulfate/	Powder	25g	Bottle	G020-25GM	Inquire	Inquire
Gentamicin Sulfate/	Powder	100g	Bottle	G020-100GM	Inquire	Inquire

Product Catalog Antibiotics Selection

Hygromycin B

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
<p>Hygromycin B/ CAS Number: 31282-04-9; Formula: C₂₀H₃₇N₃O₁₃; Formula Weight: 527. 5. Store at 2° to 8°C. This item is considered a Dangerous Good for shipping purposes and cannot be shipped to a residential address. A \$25 Dangerous Good Fee may be added per box to the shipping cost of this item.</p>	Powder	1g	Bottle	H010-1GM	Inquire	Inquire
<p>Hygromycin B/ CAS Number: 31282-04-9; Formula: C₂₀H₃₇N₃O₁₃; Formula Weight: 527. 5. Store at 2° to 8°C. This item is considered a Dangerous Good for shipping purposes and cannot be shipped to a residential address. A \$25 Dangerous Good Fee may be added per box to the shipping cost of this item.</p>	Powder	5g	Bottle	H010-5GM	Inquire	Inquire
<p>Hygromycin B/ CAS Number: 31282-04-9; Formula: C₂₀H₃₇N₃O₁₃; Formula Weight: 527. 5. Store at 2° to 8°C. This item is considered a Dangerous Good for shipping purposes and cannot be shipped to a residential address. A \$25 Dangerous Good Fee may be added per box to the shipping cost of this item.</p>	Powder	10g	Bottle	H010-10GM	Inquire	Inquire

Product Catalog Antibiotics Selection

Kanamycin

NAME/DESCRIPTION	FORMAT	SIZE	UNIT	CATALOG #	LIST PRICE	YOUR PRICE
Kanamycin Sulfate/ CAS Number: 25389-94-0; Formula: C ₁₈ H ₃₆ N ₄ O ₁₁ • H ₂ SO ₄ ; Formula Weight: 582.58. Store at 0° to -20°C.	Powder	5g	Bottle	K003-5GM	Inquire	Inquire
Kanamycin Sulfate/ CAS Number: 25389-94-0; Formula: C ₁₈ H ₃₆ N ₄ O ₁₁ • H ₂ SO ₄ ; Formula Weight: 582.58. Store at 0° to -20°C.	Powder	10g	Bottle	K003-10GM	Inquire	Inquire
Kanamycin Sulfate/ CAS Number: 25389-94-0; Formula: C ₁₈ H ₃₆ N ₄ O ₁₁ • H ₂ SO ₄ ; Formula Weight: 582.58. Store at 0° to -20°C.	Powder	25g	Bottle	K003-25GM	Inquire	Inquire
Kanamycin Sulfate/ CAS Number: 25389-94-0; Formula: C ₁₈ H ₃₆ N ₄ O ₁₁ • H ₂ SO ₄ ; Formula Weight: 582.58. Store at 0° to -20°C.	Powder	100g	Bottle	K003-100GM	Inquire	Inquire