



CERTIFICATE OF ANALYSIS

Product Name: Para-AminoBenzoic Acid

Catalog Number: A033 Lot Number: 06142010
Molecular Formula: C₇H₇NO₂
Formula Weight: 137.14 CAS Number: 150-13-0

<u>Test</u>	<u>Specifications</u>	<u>Results</u>
Appearance	White Crystalline Powder	PASS
Assay (Dried Basis)	98.0 - 102.0 %	100.7 %
Identification	Pass Test	PASS
Loss on Drying	≤ 0.2 %	0.07 %
Residue on Ignition	≤ 0.1 %	0.1%
Heavy Metals	≤ 0.002 %	<0.002 %
Aniline	≤ 10 ppm	8 ppm
p-Toluidine	≤ 10 ppm	1.1 ppm
Organic Impurities	Pass Test	PASS
Residual Solvents	Pass Test	PASS
Solubility (10%, EtOH)	Clear solution	PASS

Additional Information:

Store at 15 to 30°C

Retest Date: 31-March-2018

**This product meets or exceeds the quality standards established
by Caisson Laboratories, Inc.**

This product is intended for Laboratory Use Only

Becky Kowallis 3/23/2016
Quality Control