



CERTIFICATE OF ANALYSIS

Product Name: Banana AGS Basal Medium

Catalog Number: BMP02

Lot Number: 09130028

Storage: 2 to 8°C

Manufacture Date: SEP-2013

Expiration Date: 30-SEP-2018

MADE IN THE USA

Test	Specifications	Results
Appearance	White to Yellow Powder	PASS
Solubility in Water	4.71 g/L in ≤ 30 minutes	4 minutes
Moisture Content	≤ 2%	≤ 1.5%
Particle Size	≥ 80% through 100 mesh	> 90%
pH (measured at ~22° C)	4.25 - 5.25	4.78
Osmolality	Record mM/Kg	89 mM/Kg
Plant Tissue Culture of: <i>Solanum tuberosum</i> <i>Vitis spp.</i> <i>Musa acuminata</i> (cv.) <i>Mormata</i> <i>Musa acuminata</i> (cv.) <i>Red Dacca</i>	Supports plant growth with no morphological abnormalities and/or leaf and shoot formation.	PASS

This product was produced and tested in accordance with all relevant SOP's and meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This Product is intended for Laboratory Use Only.

Prepared by:

Kurt Rieben
Quality Control

11-JUL-2017

Reviewed by:

BSE/TSE Declaration

To the present state of our knowledge we declare that this product does not contain any BSE/TSE AGENT. This declaration is given in good faith and no warranty express or implied with respect to quality and properties is made.