



## CERTIFICATE OF ANALYSIS

Product Name: Calcium Chloride, Anhydrous  
ACS Grade

Catalog Number: C021 Lot Number: 02152021

Molecular Formula:  $\text{CaCl}_2$  CAS Number: 10043-52-4

Formula Weight: 110.99

<u>Test</u>	<u>Specifications</u>	<u>Results</u>
Appearance	White Flake-like Powder	PASS
Identification	Pass Test	PASS
Purity	>96 %	98.7 %
Titration Base	< 0.006 meq/g	0.005 meq/g
Solubility (200mg/mL, H <sub>2</sub> O)	Clear, colorless solution	PASS

**Additional Information:**

Store at 15 to 30°C

Expiration Date: Feb 28, 2019

**This product meets or exceeds the quality standards established by Caisson Laboratories, Inc.**

*This product is intended for Laboratory Use Only*

Becky Kowallis 3/12/2015  
Quality Control