



CERTIFICATE OF ANALYSIS

Product Name: Myo-inositol

Catalog Number: M014
Molecular Formula: $C_6H_{12}O_6$
Formula Weight: 180.16

Lot Number: 07152004
CAS Number: 87-89-8
Source: Plant

<u>Test</u>	<u>Specifications</u>	<u>Results</u>
Appearance	White Crystalline Powder	PASS
Identity Test	Passes Test	PASS
Assay	$\geq 99 \%$	99.50 %
Loss on Drying	$\leq 1.0 \%$	0.03 %
Solubility (50mg/mL, water)	Clear, colorless solution	PASS

Additional Information

Store at 15 to 30°C

Recommended Retest: 31-July-2019

This product meets or exceeds the quality standards established by Caisson Laboratories, Inc.

This product is intended for Laboratory Use Only

Becky Kowallis
Quality Control

8/13/2015