





Comparative Formulation

Minimum Essential Medium

 Powdered Media	MEP01	MEP02	MEP03	MEP10	MEP12	MEP14	MEP15	MEP16
Amino Acids	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
D-Valine					92.0000			
Glycine				50.0000			7.5000	7.5000
L-Alanine				25.0000			8.9000	8.9000
L-Arginine, Hydrochloride	126.0000	1260.0000	126.4000	127.0000	126.0000	126.4000	126.0000	126.0000
L-Asparagine, Monohydrate				50.0000			13.2000	13.2000
L-Aspartic Acid				30.0000			13.3000	13.3000
L-Cysteine, Hydrochloride, Monohydrate				100.0000				
L-Cystine, Dihydrochloride	32.4000	310.0000	31.2000	31.0000	31.0000	31.2000	31.0000	31.0000
L-Glutamic Acid				75.0000			14.7000	14.7000
L-Glutamine	292.0000			292.0000	292.0000	292.0000	292.0000	292.0000
L-Histidine, Hydrochloride, Monohydrate	41.8800	420.0000	41.9000	42.0000	42.0000	41.9000	42.0000	42.0000
L-Isoleucine	52.0000	520.0000	52.5000	52.0000	52.0000	52.5000	52.0000	52.0000
L-Leucine	52.0000	520.0000	52.5000	52.0000	52.0000	52.5000	52.0000	52.0000
L-Lysine, Hydrochloride	72.4600	730.0000	72.5000	73.0000	73.0000	73.0000	72.5000	72.5000
L-Methionine	15.0000	150.0000	15.0000	15.0000	15.0000	15.0000	15.0000	15.0000
L-Phenylalanine	32.0000	320.0000	32.5000	32.0000	32.0000	32.5000	32.0000	32.0000
L-Proline				40.0000			11.5000	11.5000
L-Serine				25.0000			10.5000	10.5000
L-Threonine	48.0000	480.0000	47.6000	48.0000	48.0000	47.6000	48.0000	48.0000
L-Tryptophan	10.0000	100.0000	10.0000	10.0000	10.0000	10.0000	10.0000	10.0000
L-Tyrosine, Disodium, Dihydrate	54.5200	520.0000	51.9000	52.0000	52.0000	51.9000	52.0000	52.0000
L-Valine	46.0000	460.0000	46.8000	46.0000		46.8000	46.0000	46.0000

 Powdered Media	MEP01	MEP02	MEP03	MEP10	MEP12	MEP14	MEP15	MEP16
Inorganic Salts	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
Calcium Chloride, Anhydrous		2000.0000		200.0000	200.0000	200.0000	200.0000	200.0000
Magnesium Chloride, Anhydrous	93.6800							
Magnesium Sulfate, Anhydrous		980.0000		98.0000	98.0000	97.7000	98.0000	98.0000
Potassium Chloride	400.0000	4000.0000	400.0000	400.0000	400.0000	400.0000	400.0000	400.0000
Sodium Chloride	6500.0000	68000.0000	6800.0000	6800.0000	6800.0000	6800.0000	6800.0000	6800.0000
Sodium Phosphate, Monobasic, Monohydrate	1327.0000	1400.0000	1217.0000	140.0000	140.0000	140.0000	140.0000	140.0000
Other Components	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
2'-Deoxyadenosine, Monohydrate				10.0000				
2'-Deoxycytidine, Hydrochloride				11.0000				
2'-Deoxyguanosine Monohydrate				10.0000				
Adenosine				10.0000				
Cytidine, Free Base				10.0000				
D-Glucose, Anhydrous	2000.0000	10000.0000	1000.0000	1000.0000	1000.0000	1000.0000	1000.0000	1000.0000
Guanosine				10.0000				
L-Ascorbic Acid				50.0000				
Lipoic Acid (DL-Thiotic Acid)				0.2000				
Phenol Red, Sodium Salt	10.6200			10.0000	10.6000	10.0000	10.0000	
Pyruvic Acid, Sodium Salt				110.0000				
Thymidine				10.0000				
Uridine				10.0000				



Powdered Media

MEP01

MEP02

MEP03

MEP10

MEP12

MEP14

MEP15

MEP16

Vitamins

mg/L

mg/L

mg/L

mg/L

mg/L

mg/L

mg/L

mg/L

Biotin

0.1000

Choline Chloride

1.0000

10.0000

1.0000

1.0000

1.0000

1.0000

1.0000

1.0000

D-Calcium Pantothenate

1.0000

10.0000

1.0000

1.0000

1.0000

1.0000

1.0000

1.0000

Folic Acid

1.0000

10.0000

1.0000

1.0000

1.0000

1.0000

1.0000

1.0000

Myo-Inositol

2.0000

20.0000

2.0000

2.0000

2.0000

2.0000

2.0000

2.0000

Niacinamide

1.0000

10.0000

1.0000

1.0000

1.0000

1.0000

1.0000

1.0000

Pyridoxal, Hydrochloride

1.0000

10.0000

1.0000

1.0000

1.0000

1.0000

Pyridoxine, Hydrochloride

1.0000

1.0000

Riboflavin

0.1000

1.0000

0.1000

0.1000

0.1000

0.1000

0.1000

0.1000

Thiamine, Hydrochloride

1.0000

10.0000

1.0000

1.0000

1.0000


1.0000

1.0000

1.0000

Vitamin B12

1.4000

 Powdered Media	MEP17	MEP19	MEP25	MEP29	MEP30	MEP36	MEP55
Inorganic Salts	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
Calcium Chloride, Anhydrous	200.0000	200.0000	200.0000	140.0000	140.0000	200.0000	200.0000
Ferric Chloride 6H2O			0.5400				
Magnesium Chloride, Hexahydrate			183.0000				
Magnesium Sulfate, Anhydrous	98.0000	98.0000	12.2100	98.0000	98.0000	98.0000	97.7000
Potassium Chloride	400.0000	400.0000	400.0000	400.0000	400.0000	400.0000	400.0000
Potassium Phosphate, Monobasic, Anhydrous				60.0000	60.0000		
Sodium Chloride	6800.0000	6800.0000	6800.0000	8000.0000	8000.0000	6800.0000	6800.0000
Sodium Phosphate Dibasic, Anhydrous				48.0000	48.0000		
Sodium Phosphate, Monobasic, Monohydrate	140.0000	140.0000	150.0000			140.0000	140.0000
Zinc Sulfate, Heptahydrate			0.1400				
Other Components	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
D-Glucose, Anhydrous	1000.0000	1000.0000	2000.0000	1000.0000	1000.0000	1000.0000	1000.0000
L-Ascorbic Acid	50.0000						
Linoleic Acid			0.0800				
Lipoic Acid (DL-Thiolic Acid)	0.2000		0.2000				
Meat Peptone							1500.0000
Phenol Red, Sodium Salt	10.0000	10.0000		10.0000	10.0000		10.0000
Putrescine, Dihydrochloride			0.1600				
Pyruvic Acid, Sodium Salt	110.0000		110.0000				



Powdered Media

MEP17

MEP19

MEP25

MEP29

MEP30

MEP36

MEP55

Vitamins

mg/L

mg/L

mg/L


mg/L


mg/L


mg/L


mg/L

Biotin	0.1000		0.1000					
Choline Chloride	1.0000	1.0000	56.0000	1.0000	1.0000	1.0000	1.0000	1.0000
D-Calcium Pantothenate	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Folic Acid	1.0000	1.0000	2.2000	1.0000	1.0000	1.0000	1.0000	1.0000
Myo-Inositol	2.0000	2.0000	36.0000	2.0000	2.0000	2.0000	2.0000	2.0000
Niacinamide	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Pyridoxal, Hydrochloride	1.0000	1.0000		1.0000	1.0000	1.0000		
Pyridoxine, Hydrochloride			1.0000					1.0000
Riboflavin	0.1000	0.1000	0.1000	0.1000	0.1000	0.1000	0.1000	0.1000
Thiamine, Hydrochloride	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Vitamin B12	1.4000		1.3600					

 Liquid Media	MEL01	MEL02	MEL05	MEL06	MEL07	MEL08	MEL09	MEL12
Inorganic Salts	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
Calcium Chloride, Anhydrous	2000.0000	200.0000		200.0000	200.0000	200.0000	200.0000	200.0000
Magnesium Chloride, Anhydrous			93.6800					
Magnesium Sulfate, Anhydrous	980.0000	98.0000		98.0000	98.0000	98.0000	98.0000	98.0000
Potassium Chloride	4000.0000	400.0000	400.0000	400.0000	400.0000	400.0000	400.0000	400.0000
Sodium Bicarbonate		1500.0000	2200.0000	2200.0000	2200.0000	2200.0000	2200.0000	2200.0000
Sodium Chloride	68000.0000	6800.0000	6500.0000	6800.0000	6800.0000	6800.0000	6800.0000	6800.0000
Sodium Phosphate, Monobasic, Monohydrate	1400.0000	140.0000	1327.0000	140.0000	140.0000	140.0000	140.0000	140.0000
Other Components	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
2'-Deoxyadenosine, Monohydrate					10.0000	10.0000		
2'-Deoxycytidine, Hydrochloride					11.0000	11.0000		
2'-Deoxyguanosine Monohydrate					10.0000	10.0000		
Adenosine					10.0000	10.0000		
Cytidine, Free Base					10.0000	10.0000		
D-Glucose, Anhydrous	10000.0000	1000.0000	2000.0000	1000.0000	1000.0000	1000.0000	1000.0000	1000.0000
Guanosine					10.0000	10.0000		
L-Ascorbic Acid						50.0000	50.0000	
Lipoic Acid (DL-Thiolic Acid)					0.2000	0.2000	0.2000	
Phenol Red, Sodium Salt	100.0000	10.0000	10.6200	10.0000	10.0000	10.0000	10.0000	10.6000
Pyruvic Acid, Sodium Salt					110.0000	110.0000	110.0000	
Thymidine					10.0000	10.0000		
Uridine					10.0000	10.0000		

 Liquid Media	MEL01	MEL02	MEL05	MEL06	MEL07	MEL08	MEL09	MEL12
Vitamins	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
Biotin					0.1000	0.1000	0.1000	
Choline Chloride	10.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
D-Calcium Pantothenate	10.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Folic Acid	10.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Myo-Inositol	20.0000	2.0000	2.0000	2.0000	2.0000	2.0000	2.0000	2.0000
Niacinamide	10.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Pyridoxal, Hydrochloride	10.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Riboflavin	1.0000	0.1000	0.1000	0.1000	0.1000	0.1000	0.1000	0.1000
Thiamine, Hydrochloride	10.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Vitamin B12					1.4000	1.4000	1.4000	

 Liquid Media	MEL13	MEL19	MEL20	MEL23	MEL25	MEL28	MEL32	MEL33
Inorganic Salts	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
Calcium Chloride, Anhydrous	200.0000	200.0000	200.0000	200.0000	200.0000	140.0000	140.0000	140.0000
Ferric Chloride 6H2O			0.5400		0.5400			
Magnesium Chloride, Hexahydrate			183.0000		183.0000			
Magnesium Sulfate, Anhydrous	98.0000	98.0000	12.2100	98.0000	12.2100	98.0000	98.0000	98.0000
Potassium Chloride	400.0000	400.0000	400.0000	400.0000	400.0000	400.0000	400.0000	400.0000
Potassium Phosphate, Monobasic, Anhydrous						60.0000	60.0000	60.0000
Sodium Bicarbonate	2200.0000	2200.0000	2200.0000	2200.0000	2200.0000	350.0000	350.0000	350.0000
Sodium Chloride	6800.0000	6800.0000	6800.0000	6350.0000	6800.0000	8000.0000	6800.0000	8000.0000
Sodium Phosphate Dibasic, Anhydrous						48.0000	48.0000	48.0000
Sodium Phosphate, Monobasic, Monohydrate	140.0000	140.0000	150.0000	140.0000	150.0000			
Zinc Sulfate, Heptahydrate			0.1400		0.1400			
Other Components	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
D-Glucose, Anhydrous	1000.0000	1000.0000	2000.0000	1000.0000	2000.0000	1000.0000	1000.0000	1000.0000
HEPES				5958.0000			5958.0000	
Linoleic Acid			0.0800		0.0800			
Lipoic Acid (DL-Thiotic Acid)			0.2000		0.2000			
Phenol Red, Sodium Salt	10.0000	10.0000	10.0000	10.0000		10.0000	10.0000	10.0000
Putrescine, Dihydrochloride			0.1600		0.1600			
Pyruvic Acid, Sodium Salt			110.0000		110.0000			

 Liquid Media	MEL13	MEL19	MEL20	MEL23	MEL25	MEL28	MEL32	MEL33
Vitamins	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
Biotin			0.1000		0.1000			
Choline Chloride	1.0000	1.0000	56.0000	1.0000	56.0000	1.0000	1.0000	1.0000
D-Calcium Pantothenate	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Folic Acid	1.0000	1.0000	2.2000	1.0000	2.2000	1.0000	1.0000	1.0000
Myo-Inositol	2.0000	2.0000	36.0000	2.0000	36.0000	2.0000	2.0000	2.0000
Niacinamide	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Pyridoxal, Hydrochloride	1.0000	1.0000		1.0000		1.0000	1.0000	1.0000
Pyridoxine, Hydrochloride			1.0000		1.0000			
Riboflavin	0.1000	0.1000	0.1000	0.1000	0.1000	0.1000	0.1000	0.1000
Thiamine, Hydrochloride	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000	1.0000
Vitamin B12			1.3600		1.3600			