



# Comparative Formulation

DCR

Powdered Media		DCP01
<b>Macronutrients</b>		<b>mg/L</b>
Ammonium Nitrate		400.0000
Calcium Chloride, Anhydrous		64.1400
Calcium Nitrate, Tetrahydrate		386.3100
Magnesium Sulfate, Anhydrous		180.7000
Potassium Nitrate		340.0000
Potassium Phosphate, Monobasic, Anhydrous		170.0000
<b>Micronutrients</b>		<b>mg/L</b>
Boric Acid		6.2000
Cobalt Chloride, Hexahydrate		0.0250
Cupric Sulfate, Pentahydrate		0.0250
EDTA, Disodium Salt, Dihydrate		37.3000
Ferrous Sulfate, Heptahydrate		27.8000
Manganese Sulfate, Monohydrate		22.3000
Molybdic Acid Sodium Salt, Dihydrate		0.2500
Nickel Chloride, Hexahydrate		0.0250
Potassium Iodide		0.8300
Zinc Sulfate, Heptahydrate		8.6000