




Comparative Formulation

DCR

 Powdered Media	DCP01
Macronutrients	mg/L
Ammonium Nitrate	400.0000
Calcium Chloride, Anhydrous	64.1400
Calcium Nitrate, Tetrahydrate	386.3100
Magnesium Sulfate, Anhydrous	180.7000
Potassium Nitrate	340.0000
Potassium Phosphate, Monobasic, Anhydrous	170.0000
Micronutrients	mg/L
Boric Acid	6.2000
Cobalt Chloride, Hexahydrate	0.0250
Cupric Sulfate, Pentahydrate	0.0250
EDTA, Disodium Salt, Dihydrate	37.3000
Ferrous Sulfate, Heptahydrate	27.8000
Manganese Sulfate, Monohydrate	22.3000
Molybdcic Acid Sodium Salt, Dihydrate	0.2500
Nickel Chloride, Hexahydrate	0.0250
Potassium Iodide	0.8300
Zinc Sulfate, Heptahydrate	8.6000