




Comparative Formulation

Moss

 Powdered Media	PMP01	PMP02
Macronutrients	mg/L	mg/L
Ammonium Tartrate, Dibasic	500.8880	
Calcium Nitrate, Tetrahydrate	802.9100	802.9100
Magnesium Sulfate, Anhydrous	120.3700	120.3700
Potassium Phosphate, Monobasic, Anhydrous	250.4100	250.4100
Micronutrients	mg/L	mg/L
Boric Acid	0.6140	0.6140
Cobalt Chloride, Hexahydrate	0.0547	0.0547
Cupric Sulfate, Pentahydrate	0.0549	0.0549
EDTA, Iron Sodium Salt	18.9500	18.9500
Molybdic Acid Sodium Salt, Dihydrate	0.0242	0.0242
Potassium Iodide	0.0279	0.0279
Zinc Sulfate, Heptahydrate	0.0546	0.0546