




# Comparative Formulation

LB (Miller)

 Powdered Media	LBP01	LBP02
Inorganic Salts	mg/L	mg/L
Sodium Chloride	10000.0000	10000.0000
Other Components	mg/L	mg/L
Agar, Bacteriological Type I		15000.0000
Casein Digest Peptone (Tryptone)	10000.0000	10000.0000
Yeast Extract	5000.0000	5000.0000



## Liquid Media

LBL01

Inorganic Salts

mg/L

Sodium Chloride

10000.0000

Other Components

mg/L

Casein Digest Peptone (Tryptone)

10000.0000

Yeast Extract

5000.0000