



# Comparative Formulation

## Earle's Balanced Salts

 Powdered Media	EBP01
<b>Inorganic Salts</b>	<b>mg/L</b>
Calcium Chloride, Dihydrate	264.9000
Magnesium Sulfate, Anhydrous	97.7000
Potassium Chloride	400.0000
Sodium Chloride	6800.0000
Sodium Phosphate, Monobasic, Monohydrate	140.0000
<b>Other Components</b>	<b>mg/L</b>
D-Glucose, Anhydrous	1000.0000
Phenol Red, Sodium Salt	17.0000



## Liquid Media

EBL01

### Inorganic Salts

mg/L

Calcium Chloride, Dihydrate	264.9000
Magnesium Sulfate, Anhydrous	97.7000
Potassium Chloride	400.0000
Sodium Bicarbonate	2200.0000
Sodium Chloride	6800.0000
Sodium Phosphate, Monobasic, Monohydrate	140.0000

### Other Components

mg/L

D-Glucose, Anhydrous	1000.0000
Phenol Red, Sodium Salt	17.0000