




Comparative Formulation

Phosphate-Buffered Saline

 Powdered Media	PBP01	PBP02	PBP06	PBP07
Inorganic Salts	mg/L	mg/L	mg/L	mg/L
Magnesium Chloride, Hexahydrate		100.0000		
Potassium Chloride	200.0000	200.0000		
Potassium Phosphate, Monobasic, Anhydrous	200.0000	200.0000	144.0000	210.0000
Sodium Chloride	8000.0000	8000.0000	9000.0000	9000.0000
Sodium Phosphate Dibasic, Anhydrous	1150.0000	1150.0000	421.6200	384.4200



Liquid Media

PBL01

PBL02

PBL05

PBL06

PBL07

Inorganic Salts

mg/L

mg/L

mg/L

mg/L

mg/L

Calcium Chloride, Anhydrous

100.0000

Magnesium Chloride, Hexahydrate

100.0000

Potassium Chloride

200.0000

200.0000

2000.0000

Potassium Phosphate, Monobasic, Anhydrous

200.0000

200.0000

2000.0000

144.0000

2100.0000

Sodium Chloride

8000.0000

8000.0000

80000.0000

9000.0000

90000.0000

Sodium Phosphate Dibasic, Anhydrous

1150.0000

1150.0000

11500.0000

421.6200

Sodium Phosphate, Dibasic, Heptahydrate

7260.0000