



Comparative Formulation de Fossard

 Powdered Media	DEP01	DEP02
Amino Acids	mg/L	mg/L
Glycine	3.7600	3.7600
L-Cysteine, Hydrochloride, Monohydrate	18.9200	18.9200
Macronutrients	mg/L	mg/L
Ammonium Nitrate	1601.0000	1601.0000
Calcium Chloride, Anhydrous	594.6000	594.6000
Magnesium Sulfate, Anhydrous	361.1500	361.1500
Potassium Nitrate	1011.1000	1011.1000
Sodium Phosphate, Monobasic, Monohydrate	312.0200	312.0200
Micronutrients	mg/L	mg/L
Boric Acid	18.5400	18.5400
Cobalt Chloride, Hexahydrate	0.4760	0.4760
Cupric Sulfate, Pentahydrate	0.7480	0.7480
Diammonium Molybdate	2.4800	2.4800
EDTA, Disodium Salt, Dihydrate	37.2250	37.2250
Ferrous Sulfate, Heptahydrate	27.8030	27.8030
Manganese Sulfate, Monohydrate	44.6000	44.6000
Potassium Iodide	1.6600	1.6600
Zinc Sulfate, Heptahydrate	23.0000	23.0000



Powdered Media

DEP01

DEP02

Vitamins

mg/L

mg/L

Biotin	0.2500	0.2500
Choline Chloride	1.4000	1.4000
D-Calcium Pantothenate	2.3900	2.3900
Folic Acid	0.8900	0.8900
L-Ascorbic Acid	1.7600	1.7600
Myo-Inositol	108.1000	108.1000
Nicotinic Acid	4.9300	4.9300
Pyridoxine, Hydrochloride	1.2400	1.2400
Riboflavin	3.7700	
Thiamine, Hydrochloride	13.4900	13.4900