

 Liquid Media	TRL01	TRL02	TRL03	TRL04	TRL05	TRL06
Inorganic Salts	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
EDTA, Disodium Salt, Dihydrate	372.2000			372.2000		
Potassium Chloride	400.0000			400.0000	400.0000	400.0000
Potassium Phosphate, Monobasic, Anhydrous	60.0000			60.0000	60.0000	60.0000
Sodium Bicarbonate	350.0000			350.0000	350.0000	350.0000
Sodium Chloride	8000.0000	8500.0000	8500.0000	8000.0000	8000.0000	8000.0000
Sodium Phosphate Dibasic, Anhydrous	47.5000			47.5000	47.5000	47.5000
Tetrasodium-EDTA, Anhydrous		2205.8601				220.5900
Other Components	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
D-Glucose, Anhydrous	1000.0000			1000.0000	1000.0000	1000.0000
Phenol Red, Sodium Salt	17.0000			17.0000	10.0000	10.0000
Trypsin, 1:250 from Porcine Pancreas	2500.0000	5000.0000	25000.0000	500.0000	2500.0000	500.0000