

Liquid Media	TRL01	TRL02	TRL03	TRL04	TRL05	TRL06
	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
D-Glucose, Anhydrous	1000.0000			1000.0000		
EDTA, Disodium Salt, Dihydrate	372.2000			372.2000		
Phenol Red, Sodium Salt	17.0000			17.0000		
Potassium Chloride	400.0000			400.0000		
Potassium Phosphate, Monobasic, Anhydrous	60.0000			60.0000		
Sodium Bicarbonate	350.0000			350.0000		
Sodium Chloride	8000.0000	8500.0000	8500.0000	8000.0000		
Sodium Phosphate Dibasic, Anhydrous	47.5000			47.5000		
Tetrasodium-EDTA, Anhydrous		2205.8601				220.5900
Trypsin Powder from Porcine Pancreas	2500.0000	5000.0000	25000.0000	500.0000	2500.0000	500.0000
Inorganic Salts	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
Potassium Chloride					400.0000	400.0000
Potassium Phosphate, Monobasic, Anhydrous					60.0000	60.0000
Sodium Bicarbonate					350.0000	350.0000
Sodium Chloride					8000.0000	8000.0000
Sodium Phosphate Dibasic, Anhydrous					47.5000	47.5000
Other Components	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L
D-Glucose, Anhydrous					1000.0000	1000.0000
Phenol Red, Sodium Salt					10.0000	10.0000