



# Comparative Formulation

## Quoirin & Lepoivre

 Powdered Media	QLP01	QLP03
<b>Macronutrients</b>	mg/L	mg/L
Ammonium Nitrate	400.0000	400.0000
Calcium Nitrate, Tetrahydrate	833.7700	833.7700
Magnesium Sulfate, Anhydrous	175.7900	175.7900
Potassium Nitrate	1800.0000	1800.0000
Potassium Phosphate, Monobasic, Anhydrous	270.0000	270.0000
<b>Micronutrients</b>	mg/L	mg/L
Boric Acid	6.2000	6.2000
Cobalt Chloride, Hexahydrate	0.0250	0.0250
Cupric Sulfate, Pentahydrate	0.0250	0.0250
EDTA, Disodium Salt, Dihydrate	37.3000	37.3000
Ferrous Sulfate, Heptahydrate	27.8000	27.8000
Manganese Sulfate, Monohydrate	0.7600	0.7600
Molybdcic Acid Sodium Salt, Dihydrate	0.2500	0.2500
Potassium Iodide	0.0800	0.0800
Zinc Sulfate, Heptahydrate	8.6000	8.6000
<b>Vitamins</b>	mg/L	mg/L
Myo-Inositol		100.0000
Thiamine, Hydrochloride		0.4000