



Comparative Formulation

Quoirin & Lepoivre

|  Powdered Media | QLP01 | QLP03 |
|---|-----------|-----------|
| Macronutrients | mg/L | mg/L |
| Ammonium Nitrate | 400.0000 | 400.0000 |
| Calcium Nitrate, Tetrahydrate | 833.7700 | 833.7700 |
| Magnesium Sulfate, Anhydrous | 175.7900 | 175.7900 |
| Potassium Nitrate | 1800.0000 | 1800.0000 |
| Potassium Phosphate, Monobasic, Anhydrous | 270.0000 | 270.0000 |
| Micronutrients | mg/L | mg/L |
| Boric Acid | 6.2000 | 6.2000 |
| Cobalt Chloride, Hexahydrate | 0.0250 | 0.0250 |
| Cupric Sulfate, Pentahydrate | 0.0250 | 0.0250 |
| EDTA, Disodium Salt, Dihydrate | 37.3000 | 37.3000 |
| Ferrous Sulfate, Heptahydrate | 27.8000 | 27.8000 |
| Manganese Sulfate, Monohydrate | 0.7600 | 0.7600 |
| Molybdcic Acid Sodium Salt, Dihydrate | 0.2500 | 0.2500 |
| Potassium Iodide | 0.0800 | 0.0800 |
| Zinc Sulfate, Heptahydrate | 8.6000 | 8.6000 |
| Vitamins | mg/L | mg/L |
| Myo-Inositol | | 100.0000 |
| Thiamine, Hydrochloride | | 0.4000 |