




Comparative Formulation

BT Cotton

 Powdered Media	BTP01
Macronutrients	mg/L
Calcium Chloride, Anhydrous	332.9200
Magnesium Sulfate, Anhydrous	240.7700
Potassium Nitrate	5055.0000
Potassium Phosphate, Monobasic, Anhydrous	272.0000
Micronutrients	mg/L
Boric Acid	6.1830
Cobalt Chloride, Hexahydrate	0.0240
Cupric Sulfate, Pentahydrate	0.0250
EDTA, Disodium Salt, Dihydrate	11.1700
Ferrous Sulfate, Heptahydrate	8.3400
Manganese Sulfate, Monohydrate	16.9000
Molybdic Acid Sodium Salt, Dihydrate	0.2420
Potassium Iodide	0.8300
Zinc Sulfate, Heptahydrate	8.6300
Vitamins	mg/L
Myo-Inositol	180.0000
Nicotinic Acid	0.4900
Pyridoxine, Hydrochloride	0.8200
Thiamine, Hydrochloride	1.3500