



# Comparative Formulation

BT Cotton

Powdered Media		BTP01
<b>Macronutrients</b>		<b>mg/L</b>
Calcium Chloride, Anhydrous		332.9200
Magnesium Sulfate, Anhydrous		240.7700
Potassium Nitrate		5055.0000
Potassium Phosphate, Monobasic, Anhydrous		272.0000
<b>Micronutrients</b>		<b>mg/L</b>
Boric Acid		6.1830
Cobalt Chloride, Hexahydrate		0.0240
Cupric Sulfate, Pentahydrate		0.0250
EDTA, Disodium Salt, Dihydrate		11.1700
Ferrous Sulfate, Heptahydrate		8.3400
Manganese Sulfate, Monohydrate		16.9000
Molybdc Acid Sodium Salt, Dihydrate		0.2420
Potassium Iodide		0.8300
Zinc Sulfate, Heptahydrate		8.6300
<b>Vitamins</b>		<b>mg/L</b>
Myo-Inositol		180.0000
Nicotinic Acid		0.4900
Pyridoxine, Hydrochloride		0.8200
Thiamine, Hydrochloride		1.3500