




# Comparative Formulation

White Basal

 Powdered Media	WHP01
Macronutrients	mg/L
Calcium Nitrate, Tetrahydrate	208.4700
Magnesium Sulfate, Anhydrous	351.6000
Potassium Chloride	65.0000
Potassium Nitrate	80.0000
Sodium Phosphate, Monobasic, Monohydrate	21.8500
Sodium Sulfate	200.0000
Micronutrients	mg/L
Boric Acid	1.5000
Cupric Sulfate, Pentahydrate	0.0010
Ferrous Sulfate, Heptahydrate	3.4700
Manganese Sulfate, Monohydrate	5.3100
Molybdenum Trioxide	0.0001
Potassium Iodide	0.7500
Zinc Sulfate, Heptahydrate	2.6700