



Formulations & Specifications

LB (Miller) Agar

Product Number LBP02-1LT

Luria-Bertani Medium (Miller) is a nutritionally rich bacterial growth medium. 10 g Tryptone, 5 g Yeast Extract, 10 g NaCl, 15 g Agar. Recommended usage is 40 g/L. Store at room temperature.

Components	mg/L
Agar, Bacteriological Type I	15000.0000
Casein Digest Peptone (Tryptone)	10000.0000
Sodium Chloride (NaCl)	10000.0000
Yeast Extract	5000.0000