

## Hoagland's No. 2 Basal Salt Mixture

**Product Number HOP02-50LT**

With macronutrients and micronutrients. Without ammonium phosphate. Store at 2° to 8°C.

Components	mg/L
Boric Acid ( $H_3BO_3$ )	2.8600
Calcium Nitrate, Tetrahydrate ( $Ca(NO_3)_2 \cdot 4H_2O$ )	656.4000
Cupric Sulfate, Pentahydrate ( $CuSO_4 \cdot 5H_2O$ )	0.0800
Ferric Tartrate ( $C_{12}Fe_2H_{12}O_{18}$ )	5.3200
Magnesium Sulfate, Anhydrous ( $MgSO_4$ )	240.7600
Manganese Chloride, Tetrahydrate ( $MnCl_2 \cdot 4H_2O$ )	1.8100
Molybdenum Trioxide ( $MoO_3$ )	0.0200
Potassium Nitrate ( $KNO_3$ )	606.6000
Zinc Nitrate, Hexahydrate ( $Zn(NO_3)_2 \cdot 6H_2O$ )	0.2200