



Formulations & Specifications

Trypsin-EDTA Solution 0.5% (10X)

Product Number TRL02-6X100ML

Contains 0.5% trypsin (from Porcine pancreas), 5.3mM EDTA , and 8.5 g/L NaCl. No phenol red. Sterile filtered. Porcine parvovirus and mycoplasma tested. Store at 0° to -20°C.

Components	mg/L
Sodium Chloride (NaCl)	8500.0000
Tetrasodium-EDTA, Anhydrous (C ₁₀ H ₁₂ N ₂ Na ₄)	2205.8601
Trypsin, 1:250 from Porcine Pancreas	5000.0000