

## Earle's Balanced Salts

**Product Number EBL01-6X500ML**

With calcium, magnesium, and phenol red. Sterile filtered. Store at room temperature.

Components	mg/L
Calcium Chloride, Dihydrate (CaCl <sub>2</sub> · 2H <sub>2</sub> O)	264.9000
D-Glucose, Anhydrous (C <sub>6</sub> H <sub>12</sub> O <sub>6</sub> )	1000.0000
Magnesium Sulfate, Anhydrous (MgSO <sub>4</sub> )	97.7000
Phenol Red, Sodium Salt (C <sub>19</sub> H <sub>13</sub> NaO <sub>5</sub> S)	17.0000
Potassium Chloride (KCl)	400.0000
Sodium Bicarbonate (NaHCO <sub>3</sub> )	2200.0000
Sodium Chloride (NaCl)	6800.0000
Sodium Phosphate, Monobasic, Monohydrate (NaH <sub>2</sub> PO <sub>4</sub> · H <sub>2</sub> O)	140.0000