

Earle's Balanced Salts

Product Number EBL01-6X500ML

With calcium, magnesium, and phenol red. Sterile filtered. Store at room temperature.

Components	mg/L
Calcium Chloride, Dihydrate (CaCl ₂ · 2H ₂ O)	264.9000
D-Glucose, Anhydrous (C ₆ H ₁₂ O ₆)	1000.0000
Magnesium Sulfate, Anhydrous (MgSO ₄)	97.7000
Phenol Red, Sodium Salt (C ₁₉ H ₁₃ NaO ₅ S)	17.0000
Potassium Chloride (KCl)	400.0000
Sodium Bicarbonate (NaHCO ₃)	2200.0000
Sodium Chloride (NaCl)	6800.0000
Sodium Phosphate, Monobasic, Monohydrate (NaH ₂ PO ₄ · H ₂ O)	140.0000