

HEPES-TL

Product Number IVL01-500ML

Supplement with BSA (Fract V), sodium pyruvate, and gentamicin to make HEPES-TALP (Tyrode's Albumin Lactate Pyruvate). Store at 2° to 8°C.

Components	mg/L
Calcium Chloride, Anhydrous (CaCl ₂)	220.2000
HEPES (C ₈ H ₁₃ N ₂ O ₄ S)	2384.0000
Magnesium Chloride, Hexahydrate (MgCl ₂ · 6H ₂ O)	83.9000
Potassium Chloride (KCl)	238.7200
Sodium Bicarbonate (NaHCO ₃)	168.0000
Sodium Chloride (NaCl)	6669.0000
Sodium L-Lactate solution (60%) *Measured in mL (C 3H ₅ NaO ₃)	1.4160
Sodium Phosphate, Monobasic, Monohydrate (NaH ₂ PO ₄ · H ₂ O)	48.0000