

Dulbecco's Modified Eagle's Medium (DMEM)

Product Number DMP04-10X1LT

Without Glucose, Sodium Pyruvate and amino acids.

Store at 2° to 8°C.

Components	mg/L
Calcium Chloride, Anhydrous (CaCl ₂)	200.0000
Choline Chloride (C ₅ H ₁₄ ONCl)	4.0000
D-Calcium Pantothenate (C ₁₈ H ₃₂ CaN ₂ O ₁₀)	4.0000
Ferric Nitrate, Nonahydrate (Fe(NO ₃) ₃ .9H ₂ O)	0.1000
Folic Acid (C ₁₉ H ₁₉ N ₇ O ₆)	4.0000
Magnesium Sulfate, Anhydrous (MgSO ₄)	97.6700
Myo-Inositol (C ₆ H ₁₂ O ₆)	7.2000
Niacinamide (C ₆ H ₆ N ₂ O)	4.0000
Phenol Red, Sodium Salt (C ₁₉ H ₁₃ NaO ₅ S)	15.0000
Potassium Chloride (KCl)	400.0000
Pyridoxine, Hydrochloride (C ₈ H ₁₁ NO ₃ . HCl)	4.0000
Riboflavin (C ₁₇ H ₂₀ N ₄ O ₆)	0.4000
Sodium Chloride (NaCl)	6400.0000
Sodium Phosphate, Monobasic, Monohydrate (NaH ₂ PO ₄ . H ₂ O)	125.0000
Thiamine, Hydrochloride (C ₁₂ H ₁₇ ClN ₄ OS . HCl)	4.0000