

Trypsin-EDTA Solution 0.05% (1X)

Product Number TRL04-500ML

Contains 0.05% trypsin (from Porcine pancreas) and 1mM EDTA in Hank's Balanced Salts Solution without calcium and magnesium. Sterile filtered. Porcine parvovirus and mycoplasma tested. Store at 0° to -20°C.

Components	mg/L
D-Glucose, Anhydrous (C ₆ H ₁₂ O ₆)	1000.0000
EDTA, Disodium Salt, Dihydrate (C ₁₀ H ₁₄ N ₂ Na ₂ O ₈ · 2H ₂ O)	372.2000
Phenol Red, Sodium Salt (C ₁₉ H ₁₃ NaO ₅ S)	17.0000
Potassium Chloride (KCl)	400.0000
Potassium Phosphate, Monobasic, Anhydrous (KH ₂ PO ₄)	60.0000
Sodium Bicarbonate (NaHCO ₃)	350.0000
Sodium Chloride (NaCl)	8000.0000
Sodium Phosphate Dibasic, Anhydrous (Na ₂ HPO ₄)	47.5000
Trypsin, 1:250 from Porcine Pancreas	500.0000