

Gamborg's Vitamin Mixture (Powder)

Product Number GVP01-0.1LT

Contains the vitamins as described by Gamborg et al, 1968. Each liter equivalent of powder makes one liter of a 1000X vitamin stock solution. Store at 2° to 8°C.

Components	mg/L
Myo-Inositol (C ₆ H ₁₂ O ₆)	100000.000 0
Nicotinic Acid (C ₆ H ₅ NO ₂)	1000.0000
Pyridoxine, Hydrochloride (C ₈ H ₁₁ NO ₃ . HCl)	1000.0000
Thiamine, Hydrochloride (C ₁₂ H ₁₇ CIN ₄ OS . HCl)	10000.0000