

## Amphotericin B Solution

**Product Number ABL01-100ML**

Contains 250µg/ml amphotericin B and 250µg/ml sodium deoxycholate. Effective against fungi and yeasts. Sterile-filtered solution. Store at -5° to -20°C.

Components	mg/L
Amphotericin B (C <sub>47</sub> H <sub>73</sub> NO <sub>17</sub> )	250.2500
Sodium Deoxycholate (C <sub>24</sub> H <sub>39</sub> NaO <sub>4</sub> )	250.0000