

Tryptose Phosphate Broth

Product Number TBL01-100ML

For the cultivation of a wide variety of fastidious microorganisms. Cell culture tested. Store at room temperature.

Components	mg/L
Casein Digest Peptone (Tryptone)	20000.0000
D-Glucose, Anhydrous (C ₆ H ₁₂ O ₆)	2000.0000
Sodium Chloride (NaCl)	5000.0000
Sodium Phosphate Dibasic, Anhydrous (Na ₂ HPO ₄)	2500.0000