

## Tryptose Phosphate Broth

**Product Number TBL01-100ML**

For the cultivation of a wide variety of fastidious microorganisms. Cell culture tested. Store at room temperature.

Components	mg/L
Casein Digest Peptone (Tryptone)	20000.0000
D-Glucose, Anhydrous (C <sub>6</sub> H <sub>12</sub> O <sub>6</sub> )	2000.0000
Sodium Chloride (NaCl)	5000.0000
Sodium Phosphate Dibasic, Anhydrous (Na <sub>2</sub> HPO <sub>4</sub> )	2500.0000