

Hoagland's No. 2 Basal Salt Mixture

Product Number HOP01-10LT

With macronutrients and micronutrients. Store at 2° to 8°C.

Components	mg/L
Ammonium Phosphate, Monobasic (NH ₄ H ₂ PO ₄)	115.0300
Boric Acid (H ₃ BO ₃)	2.8600
Calcium Nitrate, Tetrahydrate (Ca(NO ₃) ₂ · 4H ₂ O)	656.4000
Cupric Sulfate, Pentahydrate (CuSO ₄ · 5H ₂ O)	0.0800
Ferric Tartrate (C ₁₂ Fe ₂ H ₁₂ O ₁₈)	5.3200
Magnesium Sulfate, Anhydrous (MgSO ₄)	240.7600
Manganese Chloride, Tetrahydrate (MnCl ₂ · 4H ₂ O)	1.8100
Molybdenum Trioxide (MoO ₃)	0.0200
Potassium Nitrate (KNO ₃)	606.6000
Zinc Nitrate, Hexahydrate (Zn(NO ₃) ₂ · 6H ₂ O)	0.2200