

## Hoagland's No. 2 Basal Salt Mixture

**Product Number HOP01-10LT**

With macronutrients and micronutrients. Store at 2° to 8°C.

Components	mg/L
Ammonium Phosphate, Monobasic (NH <sub>4</sub> H <sub>2</sub> PO <sub>4</sub> )	115.0300
Boric Acid (H <sub>3</sub> BO <sub>3</sub> )	2.8600
Calcium Nitrate, Tetrahydrate (Ca(NO <sub>3</sub> ) <sub>2</sub> · 4H <sub>2</sub> O)	656.4000
Cupric Sulfate, Pentahydrate (CuSO <sub>4</sub> · 5H <sub>2</sub> O)	0.0800
Ferric Tartrate (C <sub>12</sub> Fe <sub>2</sub> H <sub>12</sub> O <sub>18</sub> )	5.3200
Magnesium Sulfate, Anhydrous (MgSO <sub>4</sub> )	240.7600
Manganese Chloride, Tetrahydrate (MnCl <sub>2</sub> · 4H <sub>2</sub> O)	1.8100
Molybdenum Trioxide (MoO <sub>3</sub> )	0.0200
Potassium Nitrate (KNO <sub>3</sub> )	606.6000
Zinc Nitrate, Hexahydrate (Zn(NO <sub>3</sub> ) <sub>2</sub> · 6H <sub>2</sub> O)	0.2200