

Hoagland's No. 2 Basal Salt Mixture

Product Number HOP01-1LT

With macronutrients and micronutrients. Store at 2° to 8°C.

Components	mg/L
Ammonium Phosphate, Monobasic ($\text{NH}_4\text{H}_2\text{PO}_4$)	115.0300
Boric Acid (H_3BO_3)	2.8600
Calcium Nitrate, Tetrahydrate ($\text{Ca}(\text{NO}_3)_2 \cdot 4\text{H}_2\text{O}$)	656.4000
Cupric Sulfate, Pentahydrate ($\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$)	0.0800
Ferric Tartrate ($\text{C}_{12}\text{Fe}_2\text{H}_{12}\text{O}_{18}$)	5.3200
Magnesium Sulfate, Anhydrous (MgSO_4)	240.7600
Manganese Chloride, Tetrahydrate ($\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$)	1.8100
Molybdenum Trioxide (MoO_3)	0.0200
Potassium Nitrate (KNO_3)	606.6000
Zinc Nitrate, Hexahydrate ($\text{Zn}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$)	0.2200