

## Hank's Balanced Salts

**Product Number HBP01-10LT**

With calcium, magnesium, and phenol red. Without sodium bicarbonate. Store at room temperature.

Components	mg/L
Calcium Chloride, Anhydrous (CaCl <sub>2</sub> )	140.0000
D-Glucose, Anhydrous (C <sub>6</sub> H <sub>12</sub> O <sub>6</sub> )	1000.0000
Magnesium Sulfate, Anhydrous (MgSO <sub>4</sub> )	98.0000
Phenol Red, Sodium Salt (C <sub>19</sub> H <sub>13</sub> NaO <sub>5</sub> S)	10.0000
Potassium Chloride (KCl)	400.0000
Potassium Phosphate, Monobasic, Anhydrous (KH <sub>2</sub> PO <sub>4</sub> )	60.0000
Sodium Chloride (NaCl)	8000.0000
Sodium Phosphate Dibasic, Anhydrous (Na <sub>2</sub> HPO <sub>4</sub> )	48.0000