

Hank's Balanced Salts

Product Number HBP01-10LT

With calcium, magnesium, and phenol red. Without sodium bicarbonate. Store at room temperature.

Components	mg/L
Calcium Chloride, Anhydrous (CaCl ₂)	140.0000
D-Glucose, Anhydrous (C ₆ H ₁₂ O ₆)	1000.0000
Magnesium Sulfate, Anhydrous (MgSO ₄)	98.0000
Phenol Red, Sodium Salt (C ₁₉ H ₁₃ NaO ₅ S)	10.0000
Potassium Chloride (KCl)	400.0000
Potassium Phosphate, Monobasic, Anhydrous (KH ₂ PO ₄)	60.0000
Sodium Chloride (NaCl)	8000.0000
Sodium Phosphate Dibasic, Anhydrous (Na ₂ HPO ₄)	48.0000