

## BT Cotton Medium

**Product Number BTP01-1LT**

With macronutrients, micronutrients, and vitamins.

Components	mg/L
Boric Acid ( $H_3BO_3$ )	6.1830
Calcium Chloride, Anhydrous ( $CaCl_2$ )	332.9200
Cobalt Chloride, Hexahydrate ( $CoCl_2 \cdot 6H_2O$ )	0.0240
Cupric Sulfate, Pentahydrate ( $CuSO_4 \cdot 5H_2O$ )	0.0250
EDTA, Disodium Salt, Dihydrate ( $C_{10}H_{14}N_2Na_2O_8 \cdot 2H_2O$ )	11.1700
Ferrous Sulfate, Heptahydrate ( $FeSO_4 \cdot 7H_2O$ )	8.3400
Magnesium Sulfate, Anhydrous ( $MgSO_4$ )	240.7700
Manganese Sulfate, Monohydrate ( $MnSO_4 \cdot H_2O$ )	16.9000
Molybdc Acid Sodium Salt, Dihydrate ( $Na_2MoO_4 \cdot 2H_2O$ )	0.2420
Myo-Inositol ( $C_6H_{12}O_6$ )	180.0000
Nicotinic Acid ( $C_6H_5NO_2$ )	0.4900
Potassium Iodide (KI)	0.8300
Potassium Nitrate ( $KNO_3$ )	5055.0000
Potassium Phosphate, Monobasic, Anhydrous ( $KH_2PO_4$ )	272.0000
Pyridoxine, Hydrochloride ( $C_8H_{11}NO_3 \cdot HCl$ )	0.8200
Thiamine, Hydrochloride ( $C_{12}H_{17}ClN_4OS \cdot HCl$ )	1.3500
Zinc Sulfate, Heptahydrate ( $ZnSO_4 \cdot 7H_2O$ )	8.6300