

BT Cotton Medium

Product Number BTP01-50LT

With macronutrients, micronutrients, and vitamins.

Components	mg/L
Boric Acid (H_3BO_3)	6.1830
Calcium Chloride, Anhydrous ($CaCl_2$)	332.9200
Cobalt Chloride, Hexahydrate ($CoCl_2 \cdot 6H_2O$)	0.0240
Cupric Sulfate, Pentahydrate ($CuSO_4 \cdot 5H_2O$)	0.0250
EDTA, Disodium Salt, Dihydrate ($C_{10}H_{14}N_2Na_2O_8 \cdot 2H_2O$)	11.1700
Ferrous Sulfate, Heptahydrate ($FeSO_4 \cdot 7H_2O$)	8.3400
Magnesium Sulfate, Anhydrous ($MgSO_4$)	240.7700
Manganese Sulfate, Monohydrate ($MnSO_4 \cdot H_2O$)	16.9000
Molybdc Acid Sodium Salt, Dihydrate ($Na_2MoO_4 \cdot 2H_2O$)	0.2420
Myo-Inositol ($C_6H_{12}O_6$)	180.0000
Nicotinic Acid ($C_6H_5NO_2$)	0.4900
Potassium Iodide (KI)	0.8300
Potassium Nitrate (KNO_3)	5055.0000
Potassium Phosphate, Monobasic, Anhydrous ($KH_2P O_4$)	272.0000
Pyridoxine, Hydrochloride ($C_8H_{11}NO_3 \cdot HCl$)	0.8200
Thiamine, Hydrochloride ($C_{12}H_{17}ClN_4OS \cdot HCl$)	1.3500
Zinc Sulfate, Heptahydrate ($ZnSO_4 \cdot 7H_2O$)	8.6300