



# Formulations & Specifications

## Gamborg's Vitamin Mixture (Powder)

**Product Number** GVP01-1LT

Contains the vitamins as described by Gamborg et al, 1968. Each liter equivalent of powder makes one liter of a 1000X vitamin stock solution. Store at 2° to 8°C.

Components	mg/L
Myo-Inositol (C <sub>6</sub> H <sub>12</sub> O <sub>6</sub> )	100000.000 0
Nicotinic Acid (C <sub>6</sub> H <sub>5</sub> NO <sub>2</sub> )	1000.0000
Pyridoxine, Hydrochloride (C <sub>8</sub> H <sub>11</sub> NO <sub>3</sub> . HCl)	1000.0000
Thiamine, Hydrochloride (C <sub>12</sub> H <sub>17</sub> ClN <sub>4</sub> OS . HCl)	10000.0000